

TSK-GEL® ODS-80T_m Products

Part Numbers:	16651, 4.6mm ID x 7.5cm, 5 μ m 08148, 4.6mm ID x 15cm, 5 μ m 08149, 4.6mm ID x 25cm, 5 μ m 14002, 21.5mm ID x 30cm, 10 μ m	19004, Guard Cartridge for 4.6mm ID columns, 3.2mm ID x 1.5cm 19018, Guard Cartridge Holder for 3.2mm ID cartridges 14098, Guard Column for P/N 14002
----------------------	---	---

This sheet contains the recommended operating conditions and the specifications for TSK-GEL ODS-80T_m columns. Installation instructions and column care information are described in the enclosed Instruction Manual.

A. OPERATING CONDITIONS	
1. Shipping Solvent:	70% Methanol - 30% Water
2. Max. Flow Rate:	1.2mL/min (4.6mm ID) 12.0mL/min (21.5mm ID) When a buffer with high viscosity is used, the maximum flow rate may have to be reduced so as not to exceed the maximum pressure drop.
3. Standard Flow Rate:	0.8 - 1.0mL/min (4.6mmID) 4.0 - 6.0mL/min (21.5mm ID)
4. Max. Pressure:	10MPa (4.6mm ID x 7.5cm column) 21MPa (4.6mm ID x 15cm column) 31MPa (4.6mm ID x 25cm column) 6MPa (21.5mm ID x 30cm)
5. pH Range:	2.0 - 7.5
6. Organic Conc.:	0 - 100%
7. Temperature:	0 - 45°C. Reduce flow rate when operating below 10°C.
8. Cleaning Solvents:	(1) Acetonitrile or Methanol, or (2) gradient from 10% to 100% acetonitrile in 0.05% trifluoroacetic acid
9. Storage:	Store the column in acetonitrile or methanol when it will not be used the next day. For overnight storage flush the column with the mobile phase at 0.2mL/min.
10. Column Protection:	The use of guard columns is recommended to prolong the life of the analytical column. Guard column life depends greatly on sample cleanliness. As a general rule, guard columns should be replaced after every 30 - 40 sample injections or when the peaks become excessively wide or when the peaks show splitting.
B. SPECIFICATIONS	
The performance of TSK-GEL ODS-80T _m columns is tested under the conditions described in the Data Sheet. All columns have passed the following quality control specifications:	
1. Number of Theoretical Plates (N):	≥ 4,500 (4.6mm ID x 7.5cm) ≥ 11,000 (4.6mm ID x 15cm) > 18,000 (4.6mm ID x 25cm) ≥ 6,000 (21.5mm ID x 30cm)
2. Asymmetry Factor (AF):	0.7 - 1.8 (4.6mm ID) 0.7 - 1.9 (21.5mm ID)