

Introduction

TOYOPEARL MX-Trp-650M is a new, high dynamic binding capacity (DBC), mixed-mode resin using tryptophan as the ligand. TOYOPEARL MX-Trp-650M contains both weak cationic and hydrophobic functional groups and is useful for protein purifications even when the target is in a high conductivity feedstock. The resin's 50-100 µm particles are stable to pressures of up to 3 bar.


Product Highlights

- High DBC (approximately 90 mg-IgG/mL-resin) when using elevated conductivity buffers
- Smaller elution pools have higher target titer and less WFI volume for use in downstream steps
- Better process throughput when run at higher linear velocities
- Excellent base stability for Clean in Place (CIP) steps

Discussion:

TOYOPEARL MX-Trp-650M is engineered to have 90-100 mg/mL DBC at higher conductivity loading conditions (Figure 1).

Figure 1. DBC Comparison at 12 mS/cm Conductivity and varying pH



Good mass transfer kinetics enable the resin to maintain its DBC at faster linear velocities (Figure 2). This fast uptake capability when coupled with narrow elution peak shape (Figure 3) results in smaller and more concentrated in-process pool volumes, reducing the amount of water for injection needed and increasing process throughput downstream. Even at conductivities of 17 mS/cm, DBC and product recovery are better than competitive mixed-mode materials (Table 1).

Figure 2. DBC at higher Linear Velocities





Figure 3. Excellent Binding and Elution Kinetics



TOYOPEARL MX-Trp-650M is base stable to 0.5 mol/L NaOH (Figure 4) making multiple uses of the resin possible after Clean in Place (CIP).


The many benefits of TOYOPEARL MX-Trp-650M shown in this overview increase process throughput and decrease cost of goods for your new therapeutic molecule.

Table 1. Recovery Comparison at Conductivities of 12 and 17 mS/cm

Resin	IgG DBC 12 mS/cm	Recovery 12 mS/cm	IgG DBC 17 mS/cm	Recovery 17 mS/cm
TOYOPEARL MX-Trp-650M	95	97%	48	96%
Capto MMC	14	86%	11	85%

Resins:	TOYOPEARL MX-Trp-650M, Capto MMC
Column Size:	6 mm ID × 4 cm
Mobile phase:	Buffer (12 mS/cm): 0.05 mol/L acetate buffer (pH 4.7) + 0.1 mol/L NaCl Buffer (17 mS/cm): 0.05 mol/L acetate buffer (pH 4.7) + 0.15 mol/L NaCl
Flow Rate:	212 cm/hr
Detection:	UV @ 280 nm
Sample:	polyclonal IgG

Figure 4. CIP Using 0.5 mol/L NaOH



Ordering Information

Part#	Description	Resin Vol.
22817	TOYOPEARL MX-Trp-650M	25 mL
22818	TOYOPEARL MX-Trp-650M	100 mL
22819	TOYOPEARL MX-Trp-650M	1 L
22820	TOYOPEARL MX-Trp-650M	5 L
22824	ToyoScreen® MX-Trp-650M	1 mL X 6
22825	ToyoScreen MX-Trp-650M	5 mL X 6

Tosoh Bioscience, TOYOPEARL and ToyoScreen are registered trademarks of Tosoh Corporation.