

Separation of Simvastatin by HPLC

Simvastatin has an action that specifically inhibits hydroxy-3-methylglutaryl-CoA (HMG-CoA) reductase in the cholesterol synthesizing system in the body, and is used as a drug to treat hyperlipidemia.

Separation of a standard sample by HPLC is shown here. Also presented is an example of highly sensitive analysis produced by modifying the eluent composition and using MS detectors.

Figure 1. Chromatogram by HPLC conditions




Table 1. Conditions

Column:	TSKgel ODS-100V, 3μm, 4.6mm ID x 5cm
Mobile phase:	0.1% H ₃ PO ₄ acetonitrile = 30/70
Flow rate:	1.0mL/min
Temperature:	25°C
Injection vol.:	5 μ L
Detection:	UV@238nm
Instrument:	Agilent 1200SL series

Figure 2. Chromatogram by MS Conditions




Table 2. Conditions

Column:	TSKgel ODS-100V, 3µm, 2.0mm ID x 15cm
Eluent:	0.1% formic acid / acetonitrile = 20/80
Gradient:	0min (10%B) → 10min (60%B) → 15min (60%B)
Flow rate:	0.2mL/min
Temperature:	40°C
Injection vol.:	2µL
Instrument:	Agilent 1200SL series QTRAP® (AB SCIEX)
Ion Source:	ESI (Positive) m/z = 419.5 (M+H)



TOSOH

TOSOH BIOSCIENCE

TOSOH Bioscience LLC
3604 Horizon Drive, Suite 100
King of Prussia, PA 19406
Orders & Service: (800) 366-4875
Fax: (610) 272-3028
www.separations.us.tosohbioscience.com
email: info.tbl@tosoh.com

**TIS136
0811**

TSKgel is a registered trademark of Tosoh Corporation.

QTrap is a registered trademark of Applied Biosystems/MDS SCIEX Instruments MDS Inc.