



Non-Porous Ion-Exchange HPLC columns for High Speed and High Resolution Analysis & Purification of Biomolecules

Hiroshi TOMIZAWA*, Kazuaki MURANAKA, Shimada MUTSUMI,
Toshiaki NISHI, Hiroyuki MORIYAMA, Yoichi YASUDA

TOSOH Corporation, Bioscience Division

www.tosohbioscience.com

Paper presented at Pittcon 2008, New Orleans; session-paper 1840-6



Abstract

- New ion-exchange HPLC columns for high-throughput protein separations were developed. They include two cation-exchange columns, one containing a sulphonic acid (SP) ligand, the other a carboxymethyl (CM) ligand. A strong anion-exchange column with quaternary ammonium groups was also developed.
- In all cases, the ligand was bonded to methacrylate-based non-porous resin particles of 5, 7 or 10 micron particle size.
- The basic properties of the new columns and their application for the analysis of biomolecules such as proteins, peptides or nucleic acids were studied.
- The pressure drops of the new ion-exchange columns were examined at various flow rates in comparison with commercial TSK-GEL DEAE-NPR and SP-NPR columns packed with 2.5 micron non-porous resin particles. As expected, all new columns exhibited substantially lower back pressures compared with the commercial NPR columns.
- Even with particle sizes above 5 micron, the new ion-exchange columns were found to provide very narrow peak widths for proteins. High-throughput separation for proteins on the ion-exchange columns was readily achievable with protein analysis times of less than 60 seconds.



Properties of Non-Porous Ion Exchange columns

- Particle chemistry: methacrylate polymer
- Particle type: non-porous
- Stationary phase ligands: Q, SP and CM
- Particle sizes of 10, 7 and 5 μ m
 - 5 μ m for Q-type only
- Column hardware and dimensions
 - Stainless steel tubing and fittings
 - 3.0mm ID x 3.5cm (10 μ m) for high-throughput separation
 - 4.6mm ID x 10cm (7 & 5 μ m) for high resolution



Benefits of Non-Porous Ion Exchange Columns

- **10 μ m: SP**
 - Low pressure drop
 - For high throughput separation of proteins, peptides and nucleic acids
 - Applicable to protein separation times of less than 60 seconds
 - Assay of reaction process for protein pegylation
- **7 μ m: SP, CM, and Q**
 - High resolution of proteins, peptides and nucleic acids
 - Separation of monoclonal antibody (mAb) aggregates
 - Separation of charged mAb isomers
 - Analysis of extend of protein pegylation
- **5 μ m: Q**
 - High resolution of nucleic acids
 - Separation of DNA oligomers
 - Separation of siRNA



Characteristics of column Types Investigated

Stationary phase	Particle size (µm)	column size (mm ID x cm)	Features	Application
Q	10	3.0 x 3.5	<ul style="list-style-type: none"> High-throughput separation of proteins, peptides and Nucleic acids Low pressure drops 	<ul style="list-style-type: none"> Proteins, peptide, LMW nucleic acids
Q	7	4.6 x 10	<ul style="list-style-type: none"> High resolution analysis of proteins, peptides 	<ul style="list-style-type: none"> Aggregates of Mab, charge isomer of mAb, PEGylated proteins
Q	5	4.6 x 10	<ul style="list-style-type: none"> High resolution analysis of nucleic acids 	<ul style="list-style-type: none"> Oligo DNA, siRNA
SP	10	3.0 x 3.5	<ul style="list-style-type: none"> High-throughput separation of proteins and peptides Low pressure drops 	<ul style="list-style-type: none"> Proteins, peptide digests
SP	7	4.6 x 10	<ul style="list-style-type: none"> High resolution analysis of proteins, peptides 	<ul style="list-style-type: none"> Aggregates of Mab, charge isomer of Mab, PEGylated proteins
CM	10	3.0 x 3.5	<ul style="list-style-type: none"> High-throughput separation of proteins and peptides Low pressure drops 	<ul style="list-style-type: none"> Proteins, peptide digests
CM	7	4.6 x 10	<ul style="list-style-type: none"> High resolution analysis of proteins, peptides 	<ul style="list-style-type: none"> Aggregates of Mab, charge isomer of mAb, PEGylated proteins



Features of Non-Porous Resin columns

Benefits of typical non-porous resin columns

- High resolution
- High speed

Challenges of typical non-porous resin columns



- High backpressure due to small particle resin
- Low loading capacity

Experimental non-porous resin columns:

- ⇒ **Much higher resolution**
- ⇒ **Faster separation (within 1 min.)**
 - **Suitable for high throughput analysis**
- ⇒ **Low backpressure**
- ⇒ **High capacity**
 - **Suitable for minor protein fraction**



Comparison of Resolution: Standard Protein Analysis by Anion Exchange columns



The same gradient condition:

Gradient : B) 0 to 100% (10min.)

The same retention time for STI:

Gradient : B) 0 to 100% (10min.) for DEAE-NPR


B) 0 to 100% (6.5min.) for Q-type trial column

Column Size: 4.6mm ID x 3.5cm, Eluent: A) 20mmol/L Tris-HCl (pH8.5), B) 0.5 mol/L NaCl in A (pH8.5)


Flow rate: 1.0mL/min, Samples: 1) conalbumin 2) ovalbumin 3) soybean trypsin inhibitor



Resolution of TSK-GEL SP-NPR Strong Cation Exchange columns with Standard Protein Analysis



High throughput analysis by
Experimental 10µm SP column




High resolution analysis by Exp.
7µm SP column

Eluent : A) 20mmol/L sodium acetate buffer (pH5.0) B) 1.0M NaCl in 20mM sodium acetate buffer (pH5.0)
Gradient : described in Figure Flow-rate : described in Figure
Sample 1) alpha-chymotrypsinogen A 2) cytochrome C 3) lysozyme Each protein 1mg/mL, 2.5uL injected



High Speed Analysis on Experimental Q column




**Analysis within 1 minute is achieved
at a flow rate of 2mL/min.**

Exp. Q column, 10 μ m, 3.0 mm ID x 3.5cm
Eluent: A) 20mmol/L Tris-HCl (pH8.5)
B) 0.5mol/L NaCl in A (pH8.5)
Gradient: B) 0% to 100% (1min.)
Flow-rate: 2.0mL/min
Sample: 1) conalbumin
2) ovalbumin
3) trypsin inhibitor
Detection: UV@280nm



Pressure as a Function of Flow Rate




Column Size: 4.6mm ID x 3.5cm
Eluent: 20mmol/L Tris-HCl (pH8.5)

Comparison of pressure drops between each particle size of the experimental columns with a conventional NPR column of the same column dimension



Comparison of loading capacity




Column Size: 4.6mm ID x 3.5cm
Eluent: A) 20mmol/L Tris-HCl (pH8.5)
B) 0.5mol/L Tris-HCl in A (pH8.5)
Flow rate: 1.0mL/min
Gradient: B) 0 to 100% (10min.)
Samples: 1) conalbumin, 2) ovalbumin
Detection: UV@ 280nm

Comparison of loading capacity between each particle size of the experimental (trial) NPR columns with a conventional NPR column of the same size

Applications of 10 μ m Non-porous Resin columns



Purification of Plasmid from *E.coli* Lysate using Experimental High Throughput Anion Exchange column



Agarose Gel Electrophoresis



Non-reduced SDS-PAGE



Exp. Q column, 10 μ m, 4.6mm ID x 3.5cm

Eluent: A) 20mmol/L Tris-HCl (pH8.5) B) 1.0mol/L NaCl in A (pH8.5)

Gradient: B) 80% to 100% (5min.) Flow-rate: 2.0 mL/min Detection: UV@280nm



Purification of PCR Products using experimental Q-type High Throughput Anion Exchange column




Agarose Gel Electrophoresis


Exp. Q column, 10 μ m, 4.6mm ID x 3.5cm
Eluent: A) 20mmol/L Tris-HCl (pH8.5), B) 1.0mol/L NaCl in A (pH8.5)
Gradient: B) 0% to 100% (1.5min.) Flow-rate: 2.0 mL/min, Detection: UV@260nm



Analysis of Nucleotides using Experimental Q-type Anion Exchange columns




Exp. Q column, $7\mu\text{m}$, $4.6\text{mm ID} \times 10\text{cm}$
Eluent: A) 20mmol/L Tris-HCl (pH8.5)
B) 0.5mol/L NaCl in A (pH8.5)
Gradient: B) 0 to 100% (10min.)
Flow-rate: 1.5mL/min.
Detection: UV@260nm



Exp. Q column, $10\mu\text{m}$, $4.6\text{mm ID} \times 3.5\text{cm}$
Eluent : A) 20mmol/L Tris-HCl (pH8.5)
B) 0.5mol/L NaCl in A (pH8.5)
Gradient : B) 0 to 100% (1min.)
Flow-rate : 4.0mL/min.
Detection : UV@260nm




Comparison of Experimental Q-type columns: High Throughput vs. High Resolution



Exp. Q column, 10µm, 4.6mmID x 3.5cm
Exp. Q column, 7µm 4.6mm ID x 10cm
Eluent: A) 20mmol/L Tris-HCl (pH8.5)
B) 0.5 mol/L NaCl in A (pH8.5)
Gradient: B) 0 to 100% (various time)
Flow-rate: 4.0mL/min for 4.6mmIDx3.5cm (10µm)
1.5mL/min for 4.6mmIDx10cm (7µm)
Sample: CDP, CTP
Detection: UV@260nm




Assay of Lysozyme PEGylation using an Experimental SP-type High Throughput Cation Exchange column



PEGylation



5mg/mL lysozyme in phosphate buffer pH 6.5
Lysozyme / PEG (Mw = 5,000) = 1 / 3



Exp. SP column, 10 μ m, 4.6mm ID x 3.5cm
Eluent: A: 20mmol/L Na acetate buffer pH5.0
B: 1.0mol/L NaCl in A pH5.0
Gradient: B: 0 to 100 % 2min
Flow-rate: 4.0mL/min
Real-time Analysis at 5-minutes intervals



Assay of beta-Lactoglobulin PEGylation using an Experimental High Throughput Cation Exchange column





PEGylation

5mg/mL beta-lactoglobulin in phosphate buffer pH6.5
Lysozyme / PEG (Mw =5,000) = 1 / 5

Exp. SP column, 10 μ m, 4.6mm ID x 3.5cm
Eluent: A: 20mmol/L Na acetate buffer pH4.5
B: 0.8mol/L NaClO₄ in A pH4.5
Gradient: B: 0 to 30 % 2min linear
Flow-rate: 4.0mL/min
Real-time Analysis at 5-minutes intervals



Assay of mAb PEGylation (1) using an Experimental High Throughput Cation Exchange column




PEGylation

5mg/mL mAb in phosphate buffer pH6.5
Lysozyme / PEG (Mw =20,000) = 1 / 5


Exp. SP column, 10 μ m, 4.6mm ID x 3.5cm
Eluent: A: 20mmol/L Na acetate buffer pH5.0
B: 1.0mol/L NaClO₄ in A pH5.0
Gradient: B: 10 to 40 % 1.8min linear
Flow-rate: 4.0mL/min
Real-time Analysis at 5-minutes intervals



Assay of mAb PEGylation (2) using an Experimental High Throughput Cation Exchange column



Tetra-PEG
Tri-PEG
Di-PEG
Mono-PEG
Non-PEG




Non-reduced SDS-PAGE



Applications of 7 μ m Non-Porous Resin columns



Comparison of Anion Exchange columns for the Analysis of BSA Digest



columns: DEAE-NPR, 2.5µm, 4.6mm ID x 3.5cm
ProPac WAX, 10µm, 4.0mm ID x 25cm
Exp. Q column, 7µm, 4.6mm ID x 10cm

Eluent: A) 20mmol/L Tris-HCl (pH8.5)
B) 0.5mol/L NaCl (pH8.5)

Gradient: B) 0 to 100% (30min.)

Flow-rate: 1.5mL/min.

Sample: BSA digest (Waters MassPREP BSA Digestion Standard)


Detection: UV@280nm



TOSOH

Comparison of Anion Exchange columns for the Separation of Mouse Ascites containing IgG (1)

Non-reduced SDS-PAGE



Fraction


Bef. After

H

**Exp. Q column
4.6mm ID x 10cm**

Mouse Ascites


Purified IgG



**ProPac WAX
4.0mm ID x 25cm**

Mouse Ascites

Purified IgG



columns: Exp. Q column, 7µm, 4.6mm ID x 10cm, ProPac WAX, 4.0mm ID x 25cm

Eluent: A) 20mmol/L Tris-HCl (pH8.5) B) 0.5mol/L NaCl in A (pH8.5)



Gradient: B) 0 to 100% (10min)

Flow rate: 1.0mL/min

Detection: UV@ 280nm




Comparison of Anion Exchange columns for the Separation of Mouse Ascites containing IgG (2)



Exp. Q column, 7 μ m, 4.6mm ID x 10cm, DEAE-NPR, 4.6mm ID x 3.5cm
Eluent: A) 20mmol/L Tris-HCl (pH8.5) B) 0.5mol/L NaCl in A (pH8.5)
Gradient: B) 0 to 100% (10min) Flow rate: 1.0mL/min Detection: UV@280nm




Simultaneous Analysis of Nucleotides using an Experimental High Resolution Anion Exchange column



Exp. Q column, 7 μ m, 4.6mm ID x 10cm
Eluent: A) 20mmol/L Tris-HCl (pH8.7)
B) 0.5 mol/L NaCl in A (pH8.7)
Gradient: B) 0 to 50% 15min Linear
Flow-rate: 1.5mL/min
Detection: UV@260nm




Simultaneous analysis of Deoxynucleotides using an Experimental High Resolution Anion Exchange column



Exp. Q column, 7µm, 4.6mm ID x 10cm
Eluent: A) 20mmol/L Tris-HCl (pH8.6)
B) 0.5 mol/L NaCl in A (pH8.6)
Gradient: B) 0 to 50% 30min Linear
Flow-rate: 1.5mL/min
Detection: UV@260nm





Simultaneous analysis of cyclic-NMP using an Experimental High Resolution Anion Exchange column



Exp. Q column, 7 μ m, 4.6mm ID x 10cm
Eluent: A) 5% AcCN in 5mmol/L Tris-HCl (pH8.6)
B) 0.15 mol/L NaCl in A (pH8.6)
Gradient: B) 0 to 100% 10min Linear
Flow-rate: 1.0mL/min
Detection: UV@260nm



Separation of PEGylated Protein using an Experimental High Resolution Cation Exchange column



Non-reduced SDS-PAGE

Exp. SP column: 7 μ m, 4.6mm ID x 10cm

Eluent: A) 20mmol/L Na acetate (pH4.2)

B) 0.5mol/L NaCl in A (pH4.2)

Gradient: B) 0 to 40% (20min.)


Sample: Somavert (Pfizer)

Flow rate: 1.5mL/min.

Detection: UV@280nm



Comparison of Anion Exchange columns for the Analysis of a mAb Sample



Eluent: A) 20mmol/L Tris-HCl (pH8.5), B) 0.5mol/L NaCl in A (pH8.5)

Gradient: B) 0 to 100% (20min.)


Flow-rate: 0.77mL/min for ProPac WAX, 1.00mL/min for Exp. Q column

Sample: MAb sample

Detection: UV@280nm






Comparison of Anion Exchange columns for the Analysis of beta-Lactoglobulin Digest



columns: DEAE-NPR, 2.5 μ m, 4.6mm ID x 3.5cm
Exp. Q column, 7 μ m, 4.6mm ID x 10cm
Eluent: A) 20mmol/L Tris-HCl (pH8.5)
B) 0.5mol/L NaCl in A (pH8.5)
Gradient: B) 0 to 33% (0-10min.)
33 to 100% (10-15min.)
Flow rate: 1.5mL/min.
Detection: UV@280nm




Comparison of Anion Exchange columns for the Analysis of mAb Digested with Pepsin (1)



Non-reduced
SDS-PAGE



Comparison of Anion Exchange columns with the Analysis of mAb Digested with Pepsin (2)




Analysis of Reaction Conversion

Exp. Q column, 7 μ m, 4.6mm ID x 10cm
Eluent: A) 20mmol/L Tris-HCl (pH8.5)
B) 0.5mol/L NaCl in A (pH8.5)
Gradient: B) 0 to 100% (10min)
Flow rate: 1.0mL/min
Detection: UV@280nm

Applications of 5 μ m Non-Porous Resin columns for Nucleic Acid Separations




Separation of 1kb DNA Ladder using an Experimental Q column



columns: Exp. Q column, 5 μ m, 4.6mm ID x 10cm, DNA-NPR, 4.6mm ID x 7.5cm
Eluent: A) 20mmol/L Tris-HCl (pH8.5). B) 1.0mol/L NaCl (pH8.5) Flow rate: 0.5mL/min
Gradient: B) 50 to 75% (20min.) for DNA-NPR, B) 75 to 100% (20min.) for Q-NPR column
Sample: 1kb Ladder Detection: UV@260nm



Separation of Same Size DNA Oligomers (26mer) using an Experimental Q column



Exp. Q column, 5 μ m, 4.6mm ID 10cm
Eluent: A) 20mmol/L Tris-HCl (pH8.5)
B) 0.75mol/L NaCl in A (pH8.5)
Gradient: B) 50 to 75% (25min)
Flow-rate: 0.8mL/min
Detection: UV@260nm

Sample: 26mer Oligomers

1. 5'-TAATTAAGGACTCCGTTCTTCTATAT-3'-NH₂
2. 5'-TCTTTACTTTAGTCACAAAGCGATAA-3'-NH₂
3. 5'-GACTCCGTTCTTCTATATTTTCGAGG-3'-NH₂
4. 5'-GGACGTGCTGGGTGTCTTCTCCGTCG-3'-NH₂



Summary

Experimental Non-Porous Resin Ion Exchange columns

- Packing Materials
 - Improved core particle
 - High resolution of biomolecules over a wide range of MW
 - Improved surface chemistry
 - Higher loading capacity and stronger retention

- Various Applications
 - High throughput analysis
 - Maintains high resolution at fast flow rates
 - Low backpressure
 - Analysis times of less than 1 minute are feasible
 - High resolution analysis
 - Much higher resolution than conventional NPR columns
 - Nucleic acid analysis
 - Excellent resolution for large DNA fragments