



Characterization and Binding Capacity Studies of PEGylated Lysozyme

Jessica CHRISTEL, Werner CONZE, Achim SPRAUER,
Volker NOEDINGER, Egbert MUELLER

Tosoh Bioscience GmbH, Stuttgart, Germany

www.tosohbioscience.com

Poster presented at HPLC 2008, Baltimore; session P-2207-W



Introduction

- PEGylation of proteins is a well established method to increase in vivo clearance rate and reduce toxicity and immunogenicity of therapeutic proteins. This polymeric modification changes the biochemical and physical properties of the protein and will therefore influence the behavior during a chromatographic purification. Resin manufacturers such as TOSOH are interested to know how a covalent PEG modification influences the performance and binding capacity of a protein on prepacked HPLC columns and chromatographic bulk media.
- Lysozyme is a well known standard protein, which is generally used to determine the dynamic binding capacity of Ion Exchange Chromatography (IEC) resins; therefore we decided to use PEG-lysozyme as a model protein in our study.
- PEGylated lysozyme is produced out of methoxy-PEG-aldehyde (with a MW of 5kDa, 10kDa and 30kDa) and hen white egg lysozyme in a phosphate buffer (0.1M, pH 6.0) in the presence of sodium-cyano-borohydrid (NaCNBH_3) as a reducing agent. The product mixture is analyzed by a TSKgel G3000SW_{xl} SEC HPLC column, SDS-PAGE, IEC (TSKgel SP-5PW (20)) and subsequent MALDI-TOF.
- Dynamic binding capacity (DBC) tests on different cation exchange media were performed with purified mono-PEGylated lysozyme.



Methods

■ PEGylation of egg white lysozyme:

- 5, 10, 30kDa methoxy PEG-aldehyde
- 100mmol/L Phosphate buffer (Na_2HPO_4 , NaH_2PO_4) pH 6.0
- PEGylation by reductive alkylation
- 20mmol/L NaCNBH_3 to reduce a Schiff base
- 100mmol/L HCl to stop PEG-reaction

■ SEC-HPLC:

- Column: TSKgel G3000SW_{XL}, 5 μm , 7.8mm ID x 30cm
- HPLC-System: Shimadzu Prominence
- Flow rate: 1mL/min
- Mobile phase: 0.1mol/L Phosphate buffer; 0.1mol/L Na_2SO_4 , pH 6.7
- Detection: UV@280nm
- Injection volume: 20 μL



Methods, Continued

■ IEC-FPLC:


- Column: TSKgel SP-5PW, 20 μ m, 6.6mm ID x 22cm, 1000Å
- Flow rate: 0.85mL/min
- Mobile phase: Buffer A: 25mmol/L Phosphate buffer; 0.1mol/L Na₂SO₄, pH 6.0
Buffer B: A + 0.5mol/L NaCl
- Detection: UV@280nm
- Injection volume: 100 μ L

■ SDS-PAGE:

- Precast gradient 4-15% Tris-HCl SDS-polyacrylamide gel
- Sample buffer: Laemmli
- Running buffer: Tris/glycine/SDS
- Staining: Silver staining




Schematic Diagram of Lysozyme




Lysozyme has six lysine residues as possible PEGylation reaction sites.




Time Course of Lysozyme PEGylation Reaction (5kDa PEG) Monitored with TSKgel G3000SW_{XL} SEC Column



1 PDA Multi 1 / 280nm 4nm



1 PDA Multi 1 / 280nm 4nm




1 PDA Multi 1 / 280nm 4nm


Di-, monoPEG and unPEGylated lysozyme species can be resolved by SEC-HPLC analysis.




Time Course of Lysozyme PEGylation Reaction (10kDa PEG) Monitored with TSKgel G3000SW_{XL} SEC Column



1 PDA Multi 1 / 280nm 4nm



1 PDA Multi 1 / 280nm 4nm




1 PDA Multi 1 / 280nm 4nm


Di-, monoPEG and unPEGylated lysozyme species can be resolved by SEC-HPLC analysis.




Time Course of Lysozyme PEGylation Reaction (30kDa PEG) Monitored with TSKgel G3000SW_{XL} SEC Column



1 PDA Multi 1 / 280nm 4nm





1 PDA Multi 1 / 280nm 4nm



1 PDA Multi 1 / 280nm 4nm



Di-, monoPEG and unPEGylated lysozyme species can be resolved by SEC-HPLC analysis.

SDS-PAGE Analysis of Lysozyme PEGylation Assay





Besides diPEG-lysozyme, monoPEG-lysozyme and unPEGylated lysozyme, a 30kDa triPEG species can be detected after 24 hours.

SDS-PAGE Analysis of Lysozyme PEGylation Assay




Besides diPEG-lysozyme, monoPEG-lysozyme and unPEGylated lysozyme, a 30kDa triPEG-lysozyme species can be detected after 24 hours.

SDS-PAGE Analysis of Lysozyme PEGylation Assay



Beside Di-, monoPEG and unPEGylated lysozyme, a 30kDa triPEG species can be detected after 24 hours.


Cation Exchange Assay of Lysozyme PEGylation: 5kDa PEG



Cation exchange assay of lysozyme PEGylation on TSKgel SP-5PW (20) resin and subsequent MALDI-TOF spectra of collected peaks.

With thanks to Dr. Harald Lange and Kai Darsow, Institute of Bioprocess Engineering, University of Erlangen, for the MALDI-TOF results.


Cation Exchange Assay of Lysozyme PEGylation: 10kDa PEG



Cation exchange assay of lysozyme PEGylation on TSKgel SP-5PW (20) resin and subsequent MALDI-TOF spectra of collected peaks.

With thanks to Dr. Harald Lange and Kai Darsow, Institute of Bioprocess Engineering, University of Erlangen, for the MALDI-TOF results.

Cation Exchange Assay of Lysozyme PEGylation: 30kDa PEG



Cation exchange assay of lysozyme PEGylation on TSKgel SP-5PW (20) resin and subsequent MALDI-TOF spectra of collected peaks.

With thanks to Dr. Harald Lange and Kai Darsow, Institute of Bioprocess Engineering, University of Erlangen, for the MALDI-TOF results.



Dynamic Binding Capacity of 1-monoPEG and Native Lysozyme on Different Cationic Exchange Resins

mg/mL					
	SP-5PW	SP-650M	SP-550C	SP-550EC	GigaCap S-650M
Particle Size	20 μ m	65 μ m	100 μ m	150 μ m	75 μ m
Nominal Pore Size	1000 Å	1000 Å	500 Å	500 Å	1000 Å
Mono-PEG lysozyme	29	18	53	28	95
Lysozyme*	33 (87%)	33 (55%)	69 (77%)	51 (54%)	121 (79%)

* Figures in brackets indicate relative % of lysozyme capacity that is accessible to mono-PEG lysozyme.



Conclusions

- SEC chromatograms show different elution volumes for monoPEG-lysozyme, diPEG-lysozyme and native lysozyme respectively for 5kDa, 10kDa and 30kDa PEG. HPLC-SEC using a TSKgel G3000SW_{XL} column is a fast tool to monitor time, concentration and temperature-dependent synthesis of PEGylated lysozyme species. An SDS-PAGE control experiment verified the SEC results.
- Cation exchange resin TSKgel SP-5PW (20) separates mono-PEGylated and di-PEGylated lysozyme isoforms, which is confirmed by MALDI-TOF analysis.
- Dynamic binding capacity studies (see **TABLE 1**) of 1-monoPEG on different cation exchange resins show significant differences. Best results regarding capacity are achieved with Toyopearl GigaCap S-650M resin; as for resolution, TSKgel SP-5PW has the best performance.
- Long term objective is to optimize LC resins (IEC and HIC) regarding high binding capacity in combination with high resolution of PEG-isoforms.