



Evaluation of High Sensitive LC/MS Analysis of Basic Compounds using a Direct Injection Column and an Ion-Suppressor Device

Seiji ITO*, Fumiya NAKATA, Yuji SAWADA

TOSOH Corporation, Kanagawa, Japan

www.tosohbioscience.com



Introduction

- Reversed phase chromatography with an ion-pair reagent in the mobile phase is generally used to obtaining sufficient retention of positively charged hydrophilic compounds by HPLC. However, when using detection by mass spectrometry, ion-pair reagents often lead to (1) a reduction in sensitivity by suppressing solute ionization, and (2) contamination of the detector.
- Hydrophilic interaction chromatography (HILIC) is frequently used for the analysis of hydrophilic compounds without using ion-pairing reagent in the mobile phase. However, this method exhibits some limitations. The chemical structure of the analyte functional groups influence such chromatographic properties as retention time, peak shape and recovery, making HILIC not as widely applicable as RPC. It has also been reported that the lifetime of HILIC columns is limited when analyzing biological samples.
- To overcome these difficulties we investigated a new, high sensitivity LC/MS method for the determination of positively charged hydrophilic compounds in biological samples and compared it with the performance of hydrophilic interaction chromatography (HILIC).



Experimental

- In this method, a suppressor device which contains strong ion exchange resin, is positioned between an analytical column and the MS detector. By this device, the ion-pair reagent is removed from the mobile phase before it enters the MS detector, thus minimizing the ion suppression effect.
- As a result, contamination of the detector is avoided and background baseline noise level and sensitivity of analytes are improved.
- A newly developed reversed phase column (TSKgel BSA-ODS-100V), in which the surface of the packing material was coated with denatured bovine serum albumin, was employed for analysis of analytes in serum samples.
- Prior to analysis, serum samples were diluted and filtered, which substantially reduce sample preparation time compared to liquid-liquid extraction or solid phase extraction (SPE) methods.



Experimental, Continued

- Four kinds of basic hydrophilic compounds (nicotine, cotinine, nornicotine and anabasine) in human serum were measured.
- The separation of these compounds was performed on a TSKgel BSA-ODS-100V (2.0mm ID x 10cm, 5 μ m) analytical column.
- Eluent A and eluent B were 10mM HFBA (heptafluoro-n-butyric acid) solution and acetonitrile, respectively.
- After injecting the serum sample, proteins present in the sample were excluded from the pores of the packing material and were switched to waste. After rotating the switching valve, the solution was directed via the suppressor device to the MS detector.
- Analytes were eluted from the column by using gradient elution and then were detected on ESI ionization.

Figure 1: Schematic Diagram of the LC/MS System






Figure 2: Protocol to Renew the Ion Exchange Resin at the Suppressor Device after each Analysis








Figure 3: Function of Ion Exchange Resin



Resin: Anion exchange resin
Resin Volume: 102μL
Flow rate: 200μL/min
Detector: Conductivity detector
Solvent: 10mM HFBA in (water/acetonitrile =1/99)
 10mM HFBA in (water/acetonitrile =30/70)
 HFBA: Heptafluoro-n-butyric acid

Removal of ion-pair reagent from mobile phase in 12 minutes

Figure 4: Function of Internal Surface of Reversed Phase Packing Material



Column: TSKgel BSA-ODS-100V (2.0mm ID x 5cm)
Eluent: 10mM HFBA in (water/acetonitrile=99/1)
Flow rate: 200 μ L/min
Detector: UV@280nm
Column temp.: 40°C
Sample vol.: 2 μ L



Conditions for Figures 5 - 7

1. RPC-Suppressor Mode

Column: TSKgel BSA-ODS-100V (2.0mm ID x 5.0cm)
Eluent A: 10mM HFBA in (water/acetonitrile = 99/1) (HFBA: Heptafluoro-n-butyric acid)
Eluent B: 10mM HFBA in acetonitrile
Flow Rate: 200 μ L/min
Column Oven: 40°C
Sample Volume: 50 μ L
Linear Gradient: 0min (B 0%), 2min (B 0%), 9min (B 30%), 11min (B 30%)
Valve Control: 0min (b), 2min (a), 11min (b)

MS: QTRAP (Applied Biosystems)

Ionization: MRM mode (ESI+)

Samples: Nicotine (163.130)

Cotinine (177.80)

Nornicotine (149.130)

Anabasine (163.118)

2. HILIC Mode

Column: TSKgel Amide-80, 3 μ m (2.0mm ID x 5.0cm)
Eluent A: 10mM ammonium formate (pH 3.75)
Eluent B: acetonitrile
Linear Gradient: 0min (B 95%), 10min (B 50%), 12min (B 50%), 13min (B 95%)
Flow Rate: 200 μ L/min
Column Oven: 40°C
Sample Volume: 10 μ L



Figure 5: Effect of Injection Volume on Chromatogram






Figure 6: Comparison of Chromatograms of Compounds Under Each Condition

(Sample: 5 μ g/L)





Nicotine

Cotinine





Nornicotine

Anabasinine

With
Suppressor



Without
suppressor



HILIC
column









Figure 7: Chromatograms of Human Serum Spiked with Standards






Table 1: Linearity, Limits of Quantification and Limits of Detection

TSKgel BSA-ODS-100V				HILIC Mode	
Calibration range $\mu\text{g/L}$	R^2	Precision, RSD $5\mu\text{g/L}$, n=6	LOQ $\mu\text{g/L}$	Calibration range $\mu\text{g/L}$	LOQ $\mu\text{g/L}$
0.2-50	0.996	0.65	0.2	1.0-50	1.5
0.1-50	0.998	0.52	0.1	No Retention	
0.1-50	0.999	0.68	0.05	1.0-50	0.6
0.1-50	0.999	0.53	0.1	1.0-50	1.5



Table 2: Precision and Recovery in the Determination of Compounds Spiked in Human Serum

	Concentration, $\mu\text{g/L}$ (Spiked Human Serum)	Precision, RSD n=6	Recovery %
Nicotine	10	1.5	92.1
	1	4.1	88.4
Cotinine	5	2.1	98.3
	0.5	2.6	96.1
Nornicotine	2	1.3	96.4
	0.2	3.1	91.5
Anabasine	5	1.6	94.7
	0.5	2.4	92.1



Results

- Deproteinization using a direct injection reversed phase column
 - Protein content in human serum was measured after separation. Recovery was more than 97.6%. This means that it is possible to analyze biological samples by direct injection without performing a lengthy deproteinization step.
- Sensitivity of method
 - Removal of HFBA from the eluent improved the background baseline noise level and peak heights of analytes became higher. Sensitivity of analytes was improved 5-25 times when compared to analysis without a suppressor device.
 - Compared with HILIC mode, with the analysis using formic acid solution as the eluent, sensitivity was improved 3 times. Furthermore, when using a TSKgel BSA-ODS-100V column, a large amount of sample could be loaded. As a result, sensitivity was improved 10-fold compared with a HILIC analysis. Additionally, cotinine was poorly retained on the HILIC column.
- Measurement of human serum
 - The level of spiked standards in human serum could be measured after simple pretreatment (10-fold dilution with eluent A followed by filtration) on this system. There were no interference peaks and recovery of all compounds was in the range of 88 to 98%.



Discussion

- By installing the ion-suppressor device before the MS detector, many compounds were analyzed using the MS detector under chromatographic conditions containing non-volatile salts or ion-pair reagents in the eluent.
- Furthermore, the ion-suppressor device improved the sensitivity of the MS detector because of a lack of contamination and ion suppression by these interferences.