



Characteristics of New Ion Chromatography Columns for High Throughput Ion Analysis

Shinji SATO, Yoshimitsu TADA and Hiroyuki MORIYAMA

Separation Center, TOSOH Corporation, Yamaguchi, Japan



Introduction

- Ion chromatography (IC) is a versatile and widely accepted technique for the analysis of inorganic and organic ions. It is characterized as a highly selective and sensitive analytical method.
- In many cases IC has superseded conventional wet chemistry methods such as titration, photometry, gravimetry, and colorimetry.
- The current trend in HPLC is high throughput analysis using ultra high pressure systems and shorter columns with smaller particles.
- To introduce this trend to IC, we have developed novel analytical columns with small particles for ion analysis.
- We evaluated the fundamental characteristics of these new TSK-GEL SuperIC-HS columns and the results are reported here. We investigated the effect of flow rate on column efficiency, selectivity as a function of column temperature, and resistance to organic solvents in addition to other characteristics.
- Furthermore, quantitative performance of ion analysis was evaluated using the ion chromatograph IC-2010 (Tosoh). Finally, we confirmed the excellent performance showing the high throughput analysis of several important environmental samples



Experimental

Columns – Tosoh Corporation (Japan)

- TSKgel SuperIC-Anion HS, 4.6mm ID x 10cm
- TSKguardcolumn SuperIC-Anion HS, 4.6mm ID x 1cm
- TSKgel SuperIC-Cation HS, 4.6mm ID x 10cm
- TSKguardcolumn SuperIC-Cation HS, 4.6mm ID x 1cm

Instrumentation

Instrument: Ion Chromatograph IC-2010 (Tosoh)

Data processing: IC-2010 WorkStation (Tosoh)

Chemicals and Reagents

All ions and reagents were purchased from Wako Chemicals (Osaka).

Water for eluent and sample preparation was purified with a Milli-Q® water purification system.

Environmental samples were taken in our office and a neighborhood river.



Separation of common ions on TSK-GEL SuperIC HS columns


Table 1: Specifications of TSK-GEL SuperIC HS columns

	TSKgel SuperIC-Anion HS	TSKgel SuperIC-Cation HS
Part No.	0022766	0022768
Column size	4.6mm ID x 10cm	4.6mm ID x 10cm
Container material	PEEK	PEEK
Gel material	Hydrophilic polymer	Poly(styrene-divinylbenzene)
Particle size	3.5 μ m	3.0 μ m
Functional group	Quaternary ammonium	Carboxylic acid
Capacity	ca. 30meq/L-gel	ca. 1.0eq/L-gel
Theoretical Plates (TP)	$\geq 8,000$ (SO_4^{2-})	$\geq 4,500$ (Na^+)



Separation of common ions on TSK-GEL SuperIC HS columns

Figure 1: Separation of 7 common anions




Columns: TSKgel SuperIC-Anion HS, 4.6mm ID x 10cm
TSKguardcolumn SuperIC-A HS, 4.6mm ID x 1cm
Mobile phase: 3.8mmol/L NaHCO₃ + 3.0mmol/L Na₂CO₃
Flow rate: 1.5mL/min
Detection: Suppressed conductivity
Temperature: 40°C
Injection vol.: 30μL
Pressure: 16MPa
Suppressor gel: TSKsuppress IC-A
Samples:
1. fluoride (1mg/L) 2. chloride (1mg/L)
3. nitrite (5mg/L) 4. bromide (5mg/L)
5. nitrate (5mg/L) 6. phosphate (10mg/l)
7. sulfate (5mg/L)



Separation of common ions on TSK-GEL SuperIC HS columns

Figure 2: Separation of 6 common cations




Columns: TSKgel SuperIC-Cation HS, 4.6mm ID x 10cm
TSKguardcolumn SuperIC-C HS, 4.6mm ID x 1cm
Mobile phase: 3.0mmol/L Methanesulfonic acid
+ 0.4mmol/L 18-crown 6-ether + 0.2mmol/L L-Histidine
Flow rate: 1.2mL/min
Detection: Suppressed conductivity
Temperature: 40°C
Injection vol.: 30 μL
Pressure: 16MPa
Suppressor gel: TSKsuppress IC-C
Samples: 1. lithium (0.5mg/L) 2. sodium (2 mg/L)
3. ammonium (2 mg/L) 4. magnesium (5 mg/L)
5. calcium (5 mg/L) 6. potassium (5 mg/L)




Evaluation of quantitative performance

Figure 3: Calibration curve of each ion



Approximation: Quadratic on each anion



Approximation: NH₄⁺ - quadratic
Other cations - linear



Evaluation of quantitative performance

Table 2: Limits of detection of ions

Anions	F ⁻	Cl ⁻	NO ₂ ⁻	Br ⁻	NO ₃ ⁻	PO ₄ ³⁻	SO ₄ ²⁻
(µg/L)	0.9	1.3	3.0	3.9	3.6	8.2	3.6
Cations	Li ⁺	Na ⁺	NH ₄ ⁺	K ⁺	Mg ²⁺	Ca ²⁺	
(µg/L)	0.1	0.3	0.3	3.0	0.2	0.4	



Analysis conditions were the same as in figures 1 and 2.

LODs were calculated by S/N = 3



Basic characteristics of TSK-GEL SuperIC HS columns


Figure 4: Ion retention as a function of column temperature





Basic characteristics of TSK-GEL SuperIC HS columns


Figure 5: Van Deemter curves





Basic characteristics of TSK-GEL SuperIC HS columns


Figure 6: Effect of organic solvents on retention and efficiency of TSKgel SuperIC-Anion HS column





Application to environmental samples

**Figure 7: Analysis of environmental samples –
Tap Water**




Analysis conditions were the same as in figure 1.



Application to environmental samples

Figure 7: Analysis of environmental samples – Tap Water, cont.




Analysis conditions were the same as in figure 1.



Application to environmental samples

**Figure 8: Analysis of environmental samples –
River Water**




Analysis conditions were the same as in figure 1.



Application to environmental samples

**Figure 8: Analysis of environmental samples –
River Water, cont.**




Analysis conditions were the same as in figure 1.



Application to environmental samples

**Figure 9: Analysis of environmental samples –
Simulated Tap Water**




Analysis conditions were the same as in figure 2.



Application to environmental samples

Figure 9: Analysis of environmental samples – Simulated Tap Water, cont.



Analysis conditions were the same as in figure 2.



Conclusions

- Common inorganic anions or cations could be analyzed in less than 5 minutes using TSKgel SuperIC-Anion HS or TSKgel SuperIC-Cation HS columns at a pressure below 20MPa.
- The TSKgel SuperIC-HS columns showed high resolution and well-balanced selectivity of common ions in an isocratic elution.
- Fast separation and detection by suppressed conductivity detection resulted in a wide dynamic range and highly sensitive detection of ions.
- As the TSKgel SuperIC-Anion HS packings are made of hydrophilic porous resin, early eluting ions such as fluoride were well-separated from the water dip.
- TSKgel SuperIC-Anion HS consist of highly crosslinked packings that are resistant to eluents containing organic solvent.
- The elution of carbonate ions in samples sometimes interferes with the quantitative determination of nitrite and bromide ions. With the use of the new TSKgel SuperIC-Anion HS, the relative position of the carbonate ions could be controlled by adjusting column temperature.