

Purification of Nucleic Acids

Table of Contents

1. Introduction	1
2. Oligonucleotides	1
3. Restriction fragments	2
4. PCR	3
5. Plasmids	4
References	5



TOSOH BIOSCIENCE LLC

TOSOH

3604 Horizon Drive, Suite 100, King of Prussia, PA 19406
Phone: (484) 805-1219 FAX: (610) 272-3028
Orders and Technical Service: (800) 366-4875
Member of the TOSOH Group

Introduction

Nucleic acids have become increasingly important in molecular biology, especially with the advent of antisense and gene therapies. To help meet this growing need, Tosoh provides a wide selection of prepacked columns and bulk media for nucleic acid analysis and purification.

Oligonucleotides


The use of oligonucleotides as primers in sequencing techniques, site-specific mutagenesis, hybridization probes, and possible therapeutic agents, like antisense, has increased the need for the synthesis, analysis and purification of high quality oligonucleotides.

A. Reversed phase separation of oligonucleotides

Most synthesis protocols for oligonucleotides incorporate the use of a protective group on the 5' terminal. Typically this protective group is dimethoxytrityl (DMT), which is a hydrophobic compound. One strategy for separating DMT-on final products from DMT failures is the use of reversed phase chromatography.

The effect of gradient conditions on the separation of 12-18mer polyadenylic oligonucleotides is shown in *Figure 1*. With the TSKgel Super-ODS column, this separation can be performed in less than five minutes under the conditions listed in *Figure 1*.

Figure 1. Separation of oligonucleotides on TSKgel Super-ODS



B. Anion exchange

Because the charge density of oligonucleotides is uniform, the net negative charge generally is a function of the length of the oligonucleotide. Thus, anion exchange separation of oligonucleotides is usually based on length. Sequence failures, therefore, can be resolved from full length products.

Figure 2 demonstrates the separation of a synthetic 20mer oligonucleotide on a TSKgel DNA-NPR column. This separation was effected using a fifteen minute linear sodium chloride gradient. As shown in *Figure 3*, retention of the oligonucleotide increased linearly with increasing column temperature². Work done with oligonucleotides on the TSKgel DEAE-5PW column demonstrated the effect of salt type and temperature on the retention of various length oligonucleotides³. *Figure 4* shows the separation of oligonucleotides on the TSKgel DEAE-5PW with an ammonium formate gradient.

Figure 2. Separation of 20mer oligonucleotide on TSKgel DNA-NPR




Figure 3. Retention time vs. Temperature





Figure 4. Separation of Oligonucleotide on TSKgel DEAE-5PW with an Ammonium Formate Gradient



Backbone-modified oligonucleotides are increasingly used for antisense therapy. These novel oligos have the benefit of longer half-lives due to resistance to endogenous nucleases. One common type of backbone-modified oligonucleotides is phosphorothioates where one of the two nonbridged oxygen atoms in the phosphate linkage has been replaced by a sulfur atom. The separation of several phosphorothioates on a TSKgel DEAE-2SW column is shown in [Figure 5](#)⁴.

Restriction fragments

Double-stranded DNA and single-stranded RNA fragments are generated when larger nucleic acid strands are treated with restriction endonucleases. The restriction endonucleases cleave the longer nucleic acids at specific sites and produce a discernible pattern of fragments. This process is used to help determine the sequence of the parent strand of DNA or RNA.

A. Ion exchange

[Figure 6](#) demonstrates a comparison of an α -DNA/Hind III digest separated on either a TSKgel DNA-NPR or DEAE-NPR column⁵. Both columns are weak ion exchangers containing diethylaminoethyl groups. However, the TSKgel DNA-NPR is longer, 7.5cm vs. 3.5cm, and has a lower dead volume than the TSKgel DEAE-NPR column. These features allow the TSKgel DNA-NPR to separate the digest more effectively. With either column, the separation is performed very rapidly.

B. Size exclusion

Size exclusion chromatography (SEC) with silica or polymer packings can be used to separate DNA and RNA fragments based on their size. DNA fragments smaller than 300 bases have been separated into discrete peaks using the TSKgel G3000SW and G4000SW columns. Recovery of the fragments from these columns was greater than 90%. A plot of chain length versus elution volume for double-stranded DNA is shown in [Figure 7](#)⁶.

Figure 5. Separation of phosphorothioates using a TSKgel DEAE-2SW





Table 1. TSKgel SEC column selection guide for DNA and RNA fragments

Molecular weight range	Base length	Recommended column
Double-stranded DNA	Base pairs	
<40,000	<55	G2000SW/SW _{XL} or G3000SW/SW _{XL}
40,000-80,000	55-110	G3000SW/SW _{XL}
80,000-250,000	110-375	G4000SW/SW _{XL}
250,000-1,000,000	375-1,500	G5000PW/PW _{XL}
RNA fragments	Base length	
<60,000	<165	G2000SW/SW _{XL} or G3000SW/SW _{XL}
60,000-120,000	165-330	G3000SW/SW _{XL}
120,000 - 1,200,000	330-3,300	G4000SW/SW _{XL}
1,200,000-10,000,000	3,300-27,500	G5000PW/PW _{XL}

Table 1 lists the SEC columns, along with the recommended molecular weight and base pair ranges, for DNA and RNA fragments⁷.


Figure 6. λ DNA/Hind III Digest Separation on TSKgel DNA-NPR and DEAE-NPR



PCR

The polymerase chain reaction (PCR) is a valuable biochemical tool which can amplify a specific piece of DNA of up to 6kb. It is used for diagnosis of disease, mutation detection and site-directed mutagenesis. Liquid chromatography is a rapid and sensitive post-PCR technique. In particular, ion exchange chromatography can be used to analyze PCR generated fragments. Figure 8 shows the detection of a 130bp target derived from HIV using a nonporous TSKgel DEAE-NPR column⁸.

Figure 7. Chain Length vs. Elution Volume for Double-Stranded DNA Fragments on TSKgel G3000SW and G4000SW Columns




Plasmids

Plasmids are extrachromosomal DNA found in bacteria and yeast which carry genetic information and replicate independently of the bacterial or yeast chromosome. They are circular molecules of duplex DNA ranging in size from 1 to 200kb. Plasmids are useful for molecular cloning. With the advent of gene therapy, there is an increased need for purified plasmids.

A. Ion Exchange


Figure 9 illustrates the separation of crude pBR322 plasmid on a TSKgel DEAE-5PW column[®]. This chromatographic separation provides purified plasmid in one hour, as opposed to a conventional Cs-Cl density gradient ultracentrifugation which can take up to three days.

Figure 8. Detection of HIV-1 PCR-Amplified 130-bp Target using a TSKgel DEAE-NPR Column



Column: TSKgel DEAE-NPR, 4.6mm ID x 3.5cm
Sample: HIV-1 PCR-amplified 130bp target
Load: 20 μ L
Flow rate: 1mL/min
Buffer A: 20mM Tris-HCl with 0.25M NaCl, pH 7.7
Buffer B: 20mM Tris-HCl with 1M NaCl, pH 7.7
Detection: UV @ 260nm
Temperature: ambient

Figure 9. Detection of HIV-1 PCR-Amplified 130-bp target using a TSKgel DEAE-NPR Column



Column: TSKgel DEAE-5PW, 7.5mm ID x 7.5 cm
Sample: pBR322 plasmid
Load: 2.5mg in 1mL
Flow rate: 1mL/min
Buffer A: 25mM Tris-HCl, 1mM EDTA, pH 8.0
Buffer B: 25mM Tris-HCl, 1mM EDTA, pH8.0, with 1M NaCl
Gradient: 25-60% B in 30 minutes
Detection: UV @ 260 nm
Temperature: ambient

References

1. Moriyama, H., et al., "Rapid Separation of Nucleotides on 2mm Porous Silica Reversed Phase Packings", Poster.
2. Fisher, J., et al., "A New Anion Exchange HPLC Column for Oligonucleotide Separations, Poster presented at the San Diego Conference", San Diego, CA, November 15-18, 1995.
3. Ozaki, H., et al., J. Chromatogr., 322 (1985) 243.
4. Kanehara, H., et al., Nucleosides and Nucleotides, 15 (1996) 399.
5. Nakatani, S., et al., "Separation of DNA Restriction Fragments by Ion-Exchange Chromatography", Poster presented at the 14th International Symposium of Proteins, Peptides, and Polynucleotides, Heidelberg, Germany, November 2-4, 1994.
6. Kato, Y., et al., J. Biochem., 95 (1984) 83.
7. Kato, Y., et al., J. Chromatogr., 266 (1983) 341.
8. Katz, E.D., and Eksteen, R., Poster presented at the 10th International Symposium of Proteins, Peptides, and Polynucleotides, Wiesbaden, Germany, 29-31 October 1990.
9. Merion, M. and Warren, W., Biotechniques, 7 (1989) 60.

Additional References

Oligonucleotides

1. Schick, C., and Martin, C., Biochem., 32 (1993) 4275.
2. Kim, R., et al., Biotechniques, 18 (1995) 992.
3. Moriyama, H., and Kato, Y., J. Chromatogr., 445 (1988) 225.
4. Baba, Y., J. Chromatogr. Biomed. Appl., 618 (1993) 41.
5. Oefner, P.J., et al., J. Chromatogr., 625 (1992) 331.
6. Kato, Y., et al., J. Chromatogr., 447 (1988) 212.
7. Yamaguchi, T., and Saneyoshi, M., Nucleosides and Nucleotides, 15 (1996) 607.
8. Koizumi, M., et al., Nucleosides and Nucleotides, 15 (1996) 505.
9. Guenther, R., et al., J. Chromatogr., 444 (1988) 79.

Restriction Fragments

1. Kato, Y., et al., J. Chromatogr., 478 (1989) 264.
2. Kato, Y., et al., J. Chromatogr., 320 (1985) 440.
3. Kato, Y., et al., J. Chromatogr., 266 (1983) 385.
4. Kasai, K., J. Chromatogr., 618 (1993) 203.
5. Waterborg, J.H., and Robertson, A.J., Nucl. Acids Res., 21 (1993) 2913.

PCR Fragments

1. Katz, E.D., J. Chromatogr., 512 (1990) 433.

Plasmids

1. Onishi, Y., et al., Anal. Biochem., 210 (1993) 63.



TOSOH

TOSOH BIOSCIENCE

TOSOH Bioscience LLC
3604 Horizon Drive, Suite 100
King of Prussia, PA 19406
Orders & Service: (800) 366-4875
Fax: (610) 272-3028
www.separations.us.tosohbioscience.com
email: info.tbi@tosoh.com