

# TSKgel NH<sub>2</sub>-100 Columns

A Stable Amino-Based Column Offering Alternate Selectivity for Polar Compounds in Hydrophilic Interaction Liquid Chromatography (HILIC)

TSKgel  
PRODUCT OVERVIEW

## Introduction

Hydrophilic interaction liquid chromatography (HILIC) is used primarily to separate polar and hydrophilic compounds that are not sufficiently retained by reversed phase. Target applications for HILIC include the analysis of hydrophilic compounds such as saccharides, glycosides, oligosaccharides, peptides, and drug metabolites.

A new HILIC column, the TSKgel NH<sub>2</sub>-100, provides a different selectivity from the well-known TSKgel Amide-80 columns. These novel amino-bonded phase columns (see Figure 1) stand out by providing much improved chemical stability, a prerequisite for achieving reproducible and reliable results. Due to a high ligand density and large surface area, these columns show high retention for very polar compounds.

TSKgel NH<sub>2</sub>-100 columns are packed with 3µm silica particles containing 100Å pores. A novel bonding strategy improves chemical stability of the bonded phase (see Figure 2). First the silica is reacted with a trimethylsilane endcapping reagent at a low stoichiometric ratio, before reaction of the remaining and accessible silanol groups with trifunctional alkylaminosilane. The resulting bonded phase provides a better safeguard against hydrolysis of the underlying silica.

## Product Highlights

- Alternative Selectivity to TSKgel Amide-80 HILIC Columns
- Novel Bonding Strategy Produces High Chemical Stability
- High Retention for Very Polar Compounds
- Ideally Suited for High Efficiency HPLC and LC/MS

Figure 1.





Figure 2.




## Applications

### Pyridylaminated Sugars

Although good separations of pyridylaminated (PA) sugars (Figure 3) can be obtained by reversed phase as well as by HILIC, all components are separated on the TSKgel NH<sub>2</sub>-100 column, while two of the components coelute on the TSKgel ODS-100V column. It is believed that all acidic PA sugars were completely separated on the TSKgel NH<sub>2</sub>-100 column due to an additional ionic interaction. The fact that the elution order of the acidic PA sugars is almost completely the opposite on the reversed phase and HILIC columns, demonstrates the potential of the TSKgel NH<sub>2</sub>-100 column as a second dimension in 2D HPLC.


Figure 3.



### MTX and its Derivatives

Polyglutamate metabolites of methotrexate (MTXPG<sub>2</sub>-7) were separated on TSKgel NH<sub>2</sub>-100 and TSKgel ODS-100V narrow bore columns (Figure 4). On the TSKgel NH<sub>2</sub>-100 column, MTX polyglutamate metabolites were eluted in the order of the number of glutamate groups in their molecules. Despite the early elution of MTXPG<sub>1</sub> and MTXPG<sub>2</sub> on the TSKgel NH<sub>2</sub>-100 column, the overall separation is superior over what can be accomplished on the C18 column.

Figure 4.



### Ordering Information

Part #	Description	Matrix	Housing	ID (mm)	Length (cm)
21967	<b>TSKgel NH<sub>2</sub>-100, 3μm</b>	Silica	Stainless Steel	2	5
21968	<b>TSKgel NH<sub>2</sub>-100, 3μm</b>	Silica	Stainless Steel	2	15
21969	<b>TSKgel NH<sub>2</sub>-100, 3μm</b>	Silica	Stainless Steel	4.6	5
21970	<b>TSKgel NH<sub>2</sub>-100, 3μm</b>	Silica	Stainless Steel	4.6	15



TOSOH

## TOSOH BIOSCIENCE

TOSOH Bioscience LLC  
 3604 Horizon Drive, Suite 100  
 King of Prussia, PA 19406  
 Orders & Service: (800) 366-4875  
 Fax: (610) 272-3028  
[www.separations.us.tosohbioscience.com](http://www.separations.us.tosohbioscience.com)  
 email: [info.tbl@tosoh.com](mailto:info.tbl@tosoh.com)