



# ***Polymer Analysis: Gel Permeation and Gel Filtration Chromatography Applications***



Joe Machamer  
GPC Marketing Manager




# Outline

- Introduction
- Calibration Curves
- Extending Linear Range
- TSKgel GPC Columns
  - TSKgel PW
  - TSKgel Alpha and SuperSW
  - TSKgel H Series
  - TSKgel SuperMultipore
- Tosoh EcoSEC GPC System



# HPLC Separation Modes






# Polymer Separation Column Types

- Water-Soluble Polymer Samples
  - TSKgel PW
  - TSKgel PW<sub>XL</sub>
  - TSKgel PW<sub>XL</sub>-CP
  - TSKgel SuperMultiporePW
- Polar Organic-Soluble Polymer Samples
  - TSKgel Alpha
  - TSKgel SuperAW
- Organic-Soluble Samples
  - TSKgel H<sub>XL</sub>
  - TSKgel SuperH
  - TSKgel H<sub>HR</sub>
  - TSKgel SuperHZ
  - TSKgel SuperMultiporeHZ


# SEC Separation Mechanism

Basis-Difference in apparent molecular size with no additional interaction between column packing material and the sample






# GPC Magic Triangle





# GPC Magic Triangle






# SEC Mechanism - No Interaction

- Proper SEC requires sample, mobile phase, and column chemistry to have similar polarities
  - Interactions between sample and column chemistry induces large errors in MW determination
  - Interactions between sample and column chemistry cause polymers to appear to be smaller






# Calibration Curve






# Resolution vs. Linear Range Trade-off

- Shallow calibration curve means:
  - More resolution
  - Less linear range



Low Resolution, Large Linear Range




Higher Resolution, Smaller Linear Range




# How To Increase Resolution

- Add more columns of same type
  - Longer Run Time
- Use smaller particle size



One column



Two columns




# How To Increase Linear Range

- Use series of columns with different pores sizes
- Use columns contained a mixture of beads with different pore sizes




# How To Increase Linear Range






# Effect on Calibration Curve





# Distortion on Chromatogram







# Schematic diagrams of SEC packings

## Conventional packings/columns


Large Pore




Medium Pore





Small Pore



Connection of conventional columns



## Multipore type packings




TSKgel SuperMultipore Type





# Comparison of Sequential vs. Multipore Columns



Columns: TSKgel G4000+3000+2500+2000H<sub>XL</sub>, 7.8mm ID x 30cm x 4  
TSKgel MultiporeH<sub>XL</sub>-M x 4

Mobile phase: THF

Flow rate: 1.0mL/min

Detection: RI

Injection vol.: 100µL

Conc.: 2.0mg/mL




# TSKgel PW Series Columns

- For water-soluble polymers
- Hydrophilic polymethacrylate stationary phase
- pH stability 2 to 12
- Organic stability to 20% Methanol
- Subtypes
  - TSKgel PW
  - TSKgel PW<sub>XL</sub>
  - TSKgel PW<sub>XL</sub>-CP
  - TSKgel SuperMultiporePW



# TSKgel PW<sub>XL</sub>

- Compared to TSKgel PW
  - Smaller Particle Size (6 - 13 $\mu$ m)
  - Higher resolution
- Pore size range <200Å to > 1000Å plus mixed bed



TSK-GEL PW<sub>XL</sub> columns: A. G2500PW<sub>XL'</sub>  
B. G3000PW<sub>XL'</sub> C. G4000PW<sub>XL'</sub> D. G5000PW<sub>XL'</sub>  
E. G6000PW<sub>XL'</sub> F. GMPW<sub>XL</sub>

Mobile phase: distilled H<sub>2</sub>O

Flow Rate: 1.0mL/min


Detection: RI

Samples: polyethylene oxides (PEO) standards  
polyethylene glycols (PEG) standards



# TSKgel PW<sub>XL</sub>-CP

- For separation of cationic water-soluble polymers
- Cationic groups on the surface prevent adsorption of cationic polymer
- Elution under low salt conditions



**TSKgel G3000PW<sub>XL</sub>-CP, 7 $\mu$ m**  
**TSKgel G5000PW<sub>XL</sub>-CP, 10 $\mu$ m**  
**TSKgel G6000PW<sub>XL</sub>-CP, 13 $\mu$ m**

Mobile phase: 0.1 mol/L NaNO<sub>3</sub>

Flow Rate: 1 mL/min


Detection: RI

Temp: 25°C

Samples: polyethylene oxides (PEO) standards  
polyethylene glycols (PEG) standards



# TSKgel PW<sub>XL</sub>-CP Cationic Polymer Analysis




TSKgel G3000PW<sub>XL</sub>-CP, 7 $\mu$ m, 7.8mm ID x 30cm  
TSKgel G5000PW<sub>XL</sub>-CP, 10 $\mu$ m, 7.8mm ID x 30cm  
TSKgel G6000PW<sub>XL</sub>-CP, 13 $\mu$ m, 7.8mm ID x 30cm

Mobile phase: 0.1mol/L NaNO<sub>3</sub>  
Flow Rate: 1mL/min  
Detection: RI  
Temperature: 25°C  
Sample Load: 3g/L, 100 $\mu$ L




# TSKgel SuperMultiporePW



- Multipore pore morphology
- Extended linear range w/o chromatogram distortion
- Less hydrophobic than TSKgel PW<sub>XL</sub> columns



# TSKgel SuperMultiporePW




**TSKgel SuperOligoPW, 6.0mm ID x 15cm**  
**TSKgel SuperMultiporePW-N, 6.0mm ID x 15cm**  
**TSKgel SuperMultiporePW-M, 6.0mm ID x 15cm**  
**TSKgel SuperMultiporePW-H, 6.0mm ID x 15cm**

Mobile phase: H<sub>2</sub>O  
Flow rate: 0.60mL/min  
Detection: RI  
Temperature: 25°C  
Samples: PEO, PEG and ethylene glycol



# PVP Analysis: TSKgel PW<sub>XL</sub> vs. TSKgel SuperMultiporePW-N



Columns: A: TSKgel SuperMultiporePW-N, 6.0mm ID x 15cm x 2  
B: TSKgel G3000PW<sub>XL</sub>+G2500PW<sub>XL</sub>, 6.0mm ID x 15cm x 2  
Mobile phase: 100mmol/L NaNO<sub>3</sub>  
Flow rate: 0.60mL/min  
Detection: RI  
Temperature: 40°C  
Injection vol.: 20µL  
Samples: PVP(K-15)





# TSKgel PW-type Applications

Sample		Column Selection		Selection Criteria	
		First selection	Second selection		
Carbohydrates	polysaccharides		TSKgel GMPW <sub>XL</sub>	TSKgel G5000PW <sub>XL</sub> + G3000PW <sub>XL</sub>	large pore size linearity of calibration curve
	oligosaccharides		TSKgel G-Oligo-PW	TSKgel G2500PW <sub>XL</sub> TSKgel G2000PW	resolving power
Nucleic Acids	DNA fragments	large	TSKgel G-DNA-PW TSKgel G5000PW <sub>XL</sub>		large pore size resolving power
		medium & small	TSKgel G4000SW TSKgel G3000SW		suitable pore size resolving power
	RNA		TSKgel G4000SW TSKgel G3000SW		
	oligonucleotides		TSKgel G2500PW <sub>XL</sub>		small pore size ionic interaction
Proteins	normal size proteins		TSKgel G3000SW TSKgel G4000SW TSKgel G2000SW	TSKgel G3000PW <sub>XL</sub> TSKgel G4000PW <sub>XL</sub>	resolving power
	large proteins	low density lipoprotein	TSKgel G6000PW <sub>XL</sub> TSKgel G5000PW <sub>XL</sub>		large pore size resolving power
		gelatin	TSKgel GMPW <sub>XL</sub>	TSKgel G5000PW <sub>XL</sub> + G3000PW <sub>XL</sub>	large pore size linearity of calibration curve
Peptides	large		TSKgel G3000SW TSKgel G2000SW	TSKgel G3000PW <sub>XL</sub>	
	small		TSKgel G25000PW <sub>XL</sub>	TSKgel G2000SW	linearity of calibration curve resolving power
Virus			TSKgel G6000PW <sub>XL</sub> TSKgel G5000PW <sub>XL</sub>		large pore size resolving power
Synthetic polymers			TSKgel GMPW <sub>XL</sub>	TSKgel G5000PW <sub>XL</sub> + G3000PW <sub>XL</sub>	large pore size linearity of calibration curve low adsorption
Synthetic oligomers	nonionic and cationic		TSKgel G-Oligo-PW	TSKgel G2500PW <sub>XL</sub>	small pore size resolving power ionic interaction
	anionic		TSKgel G2500PW <sub>XL</sub>		



# TSKgel Alpha and SuperAW Series Columns

- For Polymers of Intermediate Polarity
- Higher crosslinking than TSKgel PW series
- TSKgel SuperAW Resolution Same As Alpha, but 2x faster and 2/3 less solvent consumption
- Minimal shrink and swell in polar organic solvents
  - methanol, ACN, DMSO, IPA, THF, and HFIP
- Available in 5 discrete exclusion ranges plus a mixed bed column




# TSKgel Alpha and SuperAW Applications

- Sodium chondroitin sulfate
- Sodium alginate
- Carboxymethyl cellulose
- Sodium polystyrene sulfonate
- Polyvinyl pyrrolidone
- Gum arabic
- Ethylhydroxy-ethylcellulose
- Vinyl alcohol/vinyl butyral copolymer
- Hydroxypropylcellulose
- Polymethyl vinyl ether
- Cellulose acetate
- N-isopropyl acrylamide
- Polyacrylonitrile
- Vinyl chloride/vinyl acetate copolymer
- Styrene/allyl alcohol copolymer
- Poly (p-phenylene ether sulfone)




TOSOH

# Solvent Flexibility



**Chromatogram of Sodium Chondroitin Sulfate**

Column: TSKgel SuperAWM-H  
(6.0mm I.D. × 15cm × 2)  
Eluent: 0.2mol/L sodium nitrate  
Flow rate: 0.6mL/min  
Temperature: 40°C  
Detection: Refractive index detector  
Sample load: 20µL (0.5g/L)



**Chromatogram of Poly (p-phenylene Ether Sulfone)**

Column: TSKgel SuperAWM-H  
(6.0mm I.D. × 15cm × 2)  
Eluent: DMF containing 10mmol/L LiBr  
Flow rate: 0.6mL/min  
Temperature: 40°C  
Detection: Refractive index detector  
Sample load: 20µL (0.5g/L)



# TSKgel H-Series Columns

- For GPC Analysis of Organic Soluble Polymers
- PS-DVB Matrix
- Eight pore sizes
- Range of Mixed Bed Columns
- Range of Multipore Columns




# TSKgel H<sub>XL</sub> and SuperHZ

- TSKgel H<sub>XL</sub> (7.8mm ID x 30cm)
  - Ultra low polymer adsorption
  - Solvent can be switched one time
  - Extended linear range available with one Multipore column and five mixed-bed columns
- TSKgel SuperHZ (4.6/6.0mm ID x 30cm)
  - Smaller particle and column size versions of TSKgel H<sub>XL</sub>
  - Ultra low polymer adsorption
  - Solvent can be switched one time
  - Extended linear range available with three mixed-bed columns



# Characteristics of TSKgel MultiporeH<sub>XL</sub>-M



## Characteristics of linear type packing

- 1) The range of molecular weight for measurement is wide.
- 2) Calibration curve can be similar in a linear equation.
- 3) The distortions of chromatogram coming from the point of inflection does not appear..

**More correct measurement of molecular weight**




# TSKgel H<sub>HR</sub> and SuperH

- TSKgel H<sub>HR</sub> (7.8mm ID x 30cm)
  - Solvent can be switch many times
  - Five mixed-bed columns for extended linear range
  - Stable to 140°C
  - Broad Solvent Range
- TSKgel SuperH (6.0mm ID x 15cm)
  - Smaller particle and column size versions of TSKgel H<sub>HR</sub>
  - 2x shorter run times than TSKgel H<sub>HR</sub>





# Characteristics of Individual TSKgel H<sub>HR</sub> Columns



The range of molecular weight for measurement is narrow.

Good separation at specific range of molecular weight

TSKgel G1000H<sub>HR</sub>

Separation of solvent peaks

TSKgel G2000H<sub>HR</sub>

Oligomers

TSKgel G3000H<sub>HR</sub>

Low weight polymers

TSKgel H<sub>HR</sub> series, 7.8mm ID x 30cm

Mobile phase: THF

Flow Rate: 1.0mL/min


Detection: UV@254nm

Temperature: 25°C

Samples: polystyrene standards




# Chromatograms with Various TSKgel H<sub>HR</sub> Columns



Columns: TSKgel H<sub>HR</sub>, 7.8mm ID x 30cm x 2  
Mobile phase: THF  
Flow rate: 1.0mL/min  
Detection: RI  
Injection vol.: 100µL  
Conc.: 2.0mg/mL



# Characteristics of Mixed Bed TSKgel H<sub>HR</sub> Columns



The range of molecular weight for measurement is wide.


Good for screening measurement of molecular weight for unknown polymer

TSKgel H<sub>HR</sub> series, 7.8mm ID x 30cm

Mobile phase: THF  
Flow Rate: 1.0mL/min  
Detection: UV@254nm  
Temperature: 25°C  
Samples: polystyrene standards




# Comparison of Chromatograms with Individual and Mixed Type Packing



Columns: TSKgel H<sub>HR</sub>, 7.8mm ID x 30cm x 2  
Mobile phase: THF  
Flow rate: 1.0mL/min  
Detection: RI  
Injection vol.: 100 $\mu$ L  
Concentration: 2.0mg/mL



# Comparison of Calculated Molecular Weight with Various Type Packing



Polystyrene (NIST SRM706)

Regular type

Column: TSKgel SuperH5000+4000+3000+2000

$M_n$   $1.27 \times 10^5$   $M_w$   $2.65 \times 10^5$

Mixed type

Column: TSKgel SuperHZ-M x 4

$M_n$   $1.20 \times 10^5$   $M_w$   $2.73 \times 10^5$

Linear type

Column: TSKgel SuperMultiporeHZ-M x 4

$M_n$   $1.28 \times 10^5$   $M_w$   $2.83 \times 10^5$




# TSKgel SuperMultiporeHZ

- Ultra low polymer absorption
- Multipore pore morphology for linearity with out chromatogram distortion
- Semi-micro column size
  - Run times are 50% lower than conventional columns
  - Solvent consumption is 85% less than conventional columns




# Pore characterization of TSKgel SuperMultiporeHZ-M





# Calibration curves of TSKgel SuperMultiporeHZ columns




Column: TSKgel SuperMultiporeHZ,  
4.6mm ID x 15cm  
Mobile phase: THF  
Flow rate: 0.35 mL/min  
Detection: UV@254nm  
Temperature: 25°C  
Injection. vol.: 5µL  
Samples: Std. Polystyrene






# Properties of TSKgel SuperMultiporeHZ columns

	SuperMultiporeHZ-N	SuperMultiporeHZ-M	SuperMultiporeHZ-H
Base	PS/DVB	PS/DVB	PS/DVB
Particle Size ( $\mu\text{m}$ )	3	4	6
Exclusion Limit (MW)	120,000	2,000,000	40,000,000
Mean Pore Dia. (nm)	8	14	
Separation Range (MW)	300-50,000	500-1,000,000	1,000-10,000,000
Theoretical Plates	20,000TP/15cm	16,000TP/15cm	11,000TP/15cm
Column Size	4.6mm ID x 15cm	4.6mm ID x 15cm	4.6mm ID x 15cm





# Chromatograms of various polymers on TSKgel SuperMultiporeHZ column



Column: TSKgel SuperMultiporeHZ-H,  
4.6mm ID x 15cm x 4  
Mobile phase: THF  
Flow rate: 0.35 mL/min  
Detection: RI  
Temperature: 40°C  
Load: 10 $\mu$ L, 3g/L each  
Sample:  
1. Poly isobutylene  
2. Acrylic resin #1  
3. Polystyrene (SRM706)  
4. Polyvinylchloride  
5. Polyvinylbutyral  
6. Polycarbonate  
7. Epoxy resin  
8. Acrylic resin #2



# Semi-micro vs. Conventional Columns


	Conventional	Semi-micro
Diameter (mm)	7.8	4.6-6.0
Length (cm)	30	15
Flow Rate (mL/min.)	1.0	0.35
Run Time (Arbitrary)	100	50

Semi-micro columns reduce run times by 50% and flow rate by 2/3. Throughput is doubled and solvent costs drop by 1/6.



# Resolution: Semi-micro vs. Conventional System

- Semi-micro system offers better resolution for semi-micro columns



**TSKgel SuperHZ2000, 4.6mm ID x 15cm, x 4**

Mobile phase: THF  
Flow rate: 0.35mL/min  
Detection: RI  
Temperature: 40°C  
Injection vol.: 10 $\mu$ L (0.2mg/mL)



# EcoSEC GPC System





# Save Time and Money


- Low dead volume design and system layout
  - 7.5 $\mu$ L Stroke Volume
  - 2.5 $\mu$ L RI Cell Volume
  - Use of 0.2mm and 0.4mm tubing





# Save Time and Money

- Reduce run times by 50%



Flow rate: TSKgel SuperHZ (0.35mL/min)  
TSKgel H<sub>HR</sub> ( 1.0mL/min)

Dimensions: 4.6mm ID x 15cm\*4 (TSKgel SuperHZ)  
7.8mm ID x 30cm\*4 (TSKgel H<sub>HR</sub>)



# Save Time and Money


- Reduce solvent usage by 85%

Solvent	Competitor GPC (solvent + disposal)	EcoSEC GPC System (solvent + disposal)	<u>Savings</u>
NMP (\$30/L)	\$3082	\$1312	\$1,770.00
Chloroform (\$17/L)	\$1830	\$779	\$1,051.00
DMF (\$25/L)	\$2600.50	\$1117	\$1,463.50
HFIP (\$1000/L)	\$96,493	\$41,193	\$55,300.00






# Superior Performance



EcoSEC GPC System



Other GPC manufacturer


**TSKgel SuperMultiporeHZ-M,  
4.6mm ID x 15cm**

Mobile Phase: THF  
Flow rate: 0.35mL/min  
Detection: RI  
Temperature: 40°C  
Injection vol.: 10µL  
Sample: dicyclohexyl phthalate



# Superior Performance


- Room Temperature Affect On Retention Time
  - Conventional Instrument-No Pump Temperature Control





# Superior Performance


- Room Temperature Affect On Retention Time
  - EcoSEC GPC System-With Temperature Controlled Pumps





# Superior Performance

- Temperature Controlled Pumps
- Excellent Retention Time Precision



CV value less than 0.04% a day.

CV value less than 0.03% on different days.


TSKgel SuperMultiporeHZ-M, 4.6mm ID x 15cm, × 2

Mobile phase: THF  
Flow rate: 0.35mL/min  
Temperature: 40°C  
Injection vol.: 10 $\mu$ L  
Samples: polystyrene standards  
# of injections: 10/day



# Superior Performance

- Temperature Controlled Pumps
- Excellent Molecular Weight Precision



CV value less than 0.2% or less a day.

CV value less than 0.4% on different days.


TSKgel SuperMultiporeHZ-M, 4.6mm ID x 15cm, × 2

Mobile phase: THF  
Flow rate: 0.35mL/min  
Temperature: 40°C  
Injection vol.: 10µL  
Samples: copolymer and polyester standards  
# of injections: **10/day**



# Versatility

- Column oven holds up to 8, 30cm columns





# Conclusions

- SEC requires no interaction between sample molecules and column packing
- Columns with shallow calibration curves provide best resolution, but less linear range



# Conclusions

- Linear range can be extended
  - Individual pore size columns in sequence
  - Mixed bed column
  - Multipore column particles
    - No distortion in chromatogram





# Conclusions

- TSKgel PW Columns
  - For water-soluble polymers
- TSKgel Alpha and SuperAW Columns
  - For polar-organic soluble polymers
- TSKgel H-Series Columns
  - For organic-soluble polymers
- TSKgel SuperMultiporeHZ Columns
  - Wide linear range with no chromatogram distortion
  - Semi-micro dimensions
- Tosoh EcoSEC GPC System
  - Optimized for semi-micro columns
  - Ideal GPC system



### Available solvent of SEC columns

<b>TSKgel Column Type</b>	<b>Shipment Solvent</b>	<b>Changeable Solvent</b>
SuperHZ	Tetrahydrofuran	Benzene, Toluene, Xylene, Chloroform, Dichloroethane, Dichloromethane
H <sub>XL</sub>	Dimethylformamide	None
	Cyclohexane	None
	Acetone	n-Hexane, DMF, NMP, DMAC
SuperH HR	Tetrahydrofuran	Benzene, Toluene, Xylene, Chloroform, Dichloroethane, Dichloromethane
		Dimethylformamide, DMSO, Dioxan, n-Hexane, Cyclohexane, Dodecane, NMP, Acetone
		Quinoline, MEK, ODCB, Trichlorobenzene, HFIP, Pyridine, o-Chlorophenol/Chloroform
		Carbon Tetrachloride, Ethyl Acetate, Methanol/Chloroform, Ethanol, Dimethylacetamide
		1-Chloronaphthalene, Trichloroethane
SuperAW	Water	Methanol, Ethanol, Acetonitrile, Dimethylformamide, DMSO, Tetrahydrofuran, HFIP

Notes: In the case of TSKgel SuperHZ, only one way solvent change