



A Highly Consistent and Chemically Stable Silica-based Octadecylsilane Stationary Phase for High Performance Liquid Chromatography: TSKgel Super-ODS Column for LC/MS

André Pressley, Sheila Iuliano, Kevin O'Donnell; Tosoh Bioscience LLC, Montgomeryville, PA, USA.
Hiroshi Tomizawa, Yoshimi Hashimoto; Tosoh Corporation, Nanyo, Japan



Abstract

Octadecylsilane (ODS) columns have been most widely used for reversed phase HPLC. High throughput analysis with high sensitivity is now expected in various applications such as biological, pharmaceutical, and environmental testing. The TSK-GEL Super Series columns, composed of the TSKgel Super-ODS, TSKgel Super-Octyl, and the TSKgel Super-Phenyl are high performance columns that combine superior resolution and speed for reversed phased chromatography separations. These columns are provided in three functionalities for the selectivity required for a given separation.

The silica particles used in the Super Series columns are monodispersed spherical 2 μ m beads with a 110Å pore size. Due to the fact that the silica is high purity and metal free, there is no band broadening. In addition, the high level of endcapping minimizes secondary interactions between the sample and the column packing. The narrow particle size distribution results in low back pressure. The 2 μ m porous silica packing also provides a large surface area and higher loading capacity than similar non-porous materials. For fast separation of various compounds, the TSKGEL Super Series columns are packed in dimensions of 1mm, 2mm, and 4.6mm I.D. with 5cm and 10cm lengths. The data presented will include LC/MS separations on three different types of compounds: A) Pesticides. B) Tryptic Digest. and C) Small Peptides.



In this poster

- Comparisons between traditional 5 μ m particles (ODS-80Ts column) and ultra-fast 2 μ m particles used in the Super Series Reversed phase columns (Super-ODS, Super Octyl, and Super-Phenyl).
- Super Series columns availability in micro and semi-micro analytical format (such as 1mm and 2mm ID).
- The advantages of 2 μ m Silica.
- How to capitalize on the advantages by optimizing inter/extra column components.
- Interrelationship between sample injection volume, sample mass, flow rate, and column diameter.
- Examples of LC/MS Applications



TSK-GEL Super Series RPC Columns

Column	Functionality	Dimension	Pore size(Å)
Super-ODS	C18, polymeric 8% carbon	1, 2 or 4.6mm ID 5 or 10cm L	110
Super-Octyl	C8, polymeric 5% carbon	2 or 4.6mm ID 5 or 10cm L	110
Super-Phenyl	phenyl, polymeric 3% carbon	2 or 4.6mm ID 5 or 10cm L	110



Basic Properties of 2 μ m and 5 μ m Silica

	TSK-GEL Super Series	TSK-GEL ODS-80T _s
Silica	High purity	High purity
End-capped	Yes	Yes
Particle Size	2.2 μ m	5.06 μ m
Standard Deviation	0.27 μ m	0.87 μ m
Mean Pore Diam.	112Å	82Å
Pore Volume	0.25mL/g	0.63mL/g
Specific Surf. Area	96.8m ² /g	312.8m ² /g
Exclusion Limit	20,000Da	6,000Da
Carbon Content (C%)	~8	~15




Advantages of 2 μ m Silica

- **Reduces the diffusion path** (particle size)
- **Narrows peak bandwidth**
 - increases resolution (particle size)
 - increases sensitivity (small bore)
- **Reduces solvent consumption**
 - shorter column lengths
 - decreased analysis time
 - lower solvent strengths



Column

Optimization of Particle Size and Mobile Phase




H/u Curve in Various Eluent Compositions

Column: ● TSKgel Super-ODS (4.6mm ID x 5cm)
● TSKgel ODS-80T₈ (4.6mm ID x 15cm)
Eluent: (A) 70% methanol (B) 50% acetonitrile
Flow rate: 0.25 to 2.5mL/min
Detection: UV (254nm) Temperature: 25°C
Sample: Fluorene



Column

Optimization of Column Size: Why Smaller is Better




Relationship between Eluent Composition and Pressure

Column: A. TSKgel Super-ODS (4.6mm ID x 5cm)
B. TSKgel Super-ODS (4.6mm ID x 10cm)
Eluent: 50% CH₃OH, 50% CH₃CN
Flow rate: 0.5 to 2.0mL/min
Temperature: Ambient



Column

Effect on Column Diameter on Peak Height



Column: TSKgel Super-ODS
(A) 1.0mm I.D. _ 5cm
(B) 2.0mm I.D. _ 5cm
(C) 4.6mm I.D. _ 5cm

Eluent; A: H₂O/CH₃CN(95/5)+0.1% TFA
B: H₂O/CH₃CN(50/50)+0.1% TFA

Gradient: 1-100% B(15min, Linear gradient)



Flow rate: (A) 0.1mL/min
(B) 0.4mL/min
(C) 2.12mL/min

Inj. Vol.: 0.5μL
Temp.: 40°C
Sample: β-Lactoglobulin tryptic digest



Conditions

Effect of Flow Rate on Resolution (gradient volume is constant)





Column: TSKgel Super-ODS (2.0mm I.D. x 5cm)
Eluent: A: H₂O/CH₃CN(95/5)+0.1% TFA
B: H₂O/CH₃CN(50/50)+0.1% TFA
Gradient: 0-100% B((A) 2.5min, (B) 5min, (C) 20min, Linear gradient)
Flow rate: (A) 0.44cm/min, (B) 0.22cm/min, (C) 0.055cm/min
Inj. Vol.: 2μL
Temp.: 40°C
UV: 215nm
Sample: β-Lactoglobulin tryptic digest



Conditions

Effect of Sample Mass on Resolution




Column: TSKgel Super-ODS (2.0mm I.D. x 5cm)
Eluent: A: H₂O/CH₃CN(95/5)+0.1% TFA
B: H₂O/CH₃CN(50/50)+0.1% TFA
Gradient: 0-100% B((A) 30min, (B) 15min, (C) 5min, Linear gradient)
Flow rate: 0.4mL/min
Inj. Vol.: 2µL
Temp.: 40°C
UV: 215nm
Sample: β-Lactoglobulin tryptic digest
(A) 44.5mg/mL, (B) 22.5mg/mL, (C) 2.25mg/mL



Conditions

Sample Injection Volume and Column Efficiency (Theoretical Plates)



Column: (A), (B) TSKgel Super-ODS (4.6mm I.D.x 5cm)
Eluent: 70% methanol
Flow rate: 1.0mL/min
Temperature: 25°C
Detection: UV (254nm)
Samples: Naphthalene (0.1g/L) dissolved in (A) 70% methanol, (B) dissolved in 40% methanol, (0.1 g/L)



Detector / Tubing

Optimization of Tubing and Its Effect on Column Efficiency

Injector/Column*			Column/Detector**		
Length of the tubing	Volume of tubing	HETP	Length of the tubing	Volume of tubing	HETP
(cm)	(μ L)	(μ m)	(cm)	(μ L)	(μ m)
10.00	0.79	4.66	10.00	0.79	4.66
15.00	1.19	4.70	15.00	1.19	4.70
30.00	2.36	5.23	30.00	2.36	4.74
50.00	3.93	5.51	50.00	3.93	5.35
70.00	5.50	5.89	70.00	5.50	5.54

* Distance between Injector/column
** Distance between Column/Detector

Column: TSKgel Super-ODS (4.6mm x 5cm)
Eluent: 70% methanol Flow rate: 1mL/min
Detection: UV (254nm), micro flow cell
Sample: Fluorene



Detector / Tubing

Relationship between Detector Response and Column Efficiency

Time Constant	Naphthalene Theoretical Plates TP/column (relative reduction)	Resolution Tol/Nap
50 msec	10529 (0%)	13.37
1 sec	6996 (34%)	10.37
3 sec	3420 (68%)	6.87
Eluent:	70% Methanol	
Samples:	Toluene (Tol), Naphthalene (Nap)	



Detector / Tubing

Optimization of Detector Cell Volume and Its Effect on Column Efficiency

Cell Volume (μL)	Type	Column theoretical plates (relative reduction in theoretical plates) TP/5cm column
2	micro flow cell	10769 (0%)
10	low dead volume type	10150 (6%)
10	standard flow cell	3104 (71%)


Low dead volume means heat sink was removed.
Standard flow cell means heat sink was included.

Eluent: 70% Methanol
Sample: Fluorene



Detector / Tubing

Effect of Detector Time Constant on Theoretical Plates




Column: TSKgel Super-ODS (4.6mm ID x 5cm)
Eluent: 70% Methanol
Flow rate: 1.0mL/min
Sample: Fluorene Temperature: 25°C
Detector: UV (254nm), micro flow cell
Time constant: (A) 50 msec, (B) 1 sec, (C) 3 sec



Applications

LC/MS Chromatogram (TIC) of β -lactoglobulin Tryptic Digest





Column: TSKgel Super-ODS (2.0mm I.D. x 5cm)
Eluent: A: H₂O/CH₃CN(95/5)+0.1% TFA
B: H₂O/CH₃CN(50/50)+0.1% TFA
Gradient: 0-100% B(15 min. Linear gradient)
Flow rate: 0.4mL/min
Inj. Vol.: 2 μ L
Temp.: 40°C
UV: 215nm
Sample: β -Lactoglobulin tryptic digest
Ion mode: ESI+
Ion species: Normal Ion [MF-Liner]
TIC Range: m/z 50 to 2600



Applications

LC/MS Chromatogram Pesticides Extracted from Orange



Column: TSKgel Super-ODS (2.0mmI.D.x10cm)
Flow rate: 0.4mL/min
Inj. Vol.: 2µL
Eluent: (A) H₂O+0.1% TFA
(B) CH₃CN+0.1% TFA
Temp.: 40°C
UV: 254nm
Gradient: 5% B-100% B(12min, Linear gradient)
Ion mode: ESI(-)
Sample: Orange extracts spiked std. pesticides 1. Azimsulfron 2. Flazasulfron 3. Halosulfuron-methyl
Monitored ions: 406, 423, 433, 435, 835, 869, 889
TIC Range: m/z 100 to 900



Conclusions

- The advantages of the Super Series Columns (2 μ m particles) over the more traditional 5 μ m particles are:
 - Increased efficiency
 - Reduction in analysis time
 - Higher resolution
- The Super Series columns are available in:
 - 1, 2, or 4.6mm ID
 - 5 or 10cm Length
- The advantages of 2 μ m silica:
 - Reduction of diffusion path caused by a narrowing of particle interstices.
 - Increased resolution due to small particle size and narrow pore size distribution.
- Optimization of the inter/extra column components (such as void volume, detector response, and sample injection volume) significantly increases column performance.
- These ultra efficient columns are coupled with the specificity of Mass Spec, which results in superior analytical power.