



Determination of Hydrophilic Drugs and their Metabolites in Blood using a TSKgel HILIC column

Hiroyuki MORIYAMA, Hiroyuki YAMASAKI, Michiko SAKATA,
and Yasutoshi KAWAI

Separation Center, TOSOH Corporation, Yamaguchi, Japan



Introduction

- Reversed Phase Chromatography (RPC) and Hydrophilic Interaction Chromatography (HILIC) columns are widely utilized for separations of pharmaceutical products, which have hydrophobic and hydrophilic properties. HILIC columns possessing various hydrophilic ligands, such as hydroxyl, carbamoyl, amino and ionic groups, strongly retain hydrophilic and ionic compounds so that different selectivity can be expected in comparison with ODS columns.
- Therapeutic drug monitoring (TDM) by HPLC is required for several drug classes as the effective dosage range, varying between insufficient activity and toxic levels, is patient dependent. Although an ODS column is adopted for separating drugs in blood samples, hydrophilic compounds show poor retention times on an RPC column. We investigated the possibility of the separation of hydrophilic drugs and its metabolites in blood on a new HILIC column.
- In this poster the solvent composition and organic solvent effects on elution behavior for a new HILIC column from Tosoh Corporation, a TSK-GEL NH₂-100, 3μm column, is demonstrated. The separation of hydrophilic drugs and its metabolites in serum by on-line and off-line deproteinization procedures are also reported.



Experimental

Columns – Tosoh Corporation (Japan)

- TSKgel NH₂-100, 3μm, 4.6mm ID x 15cm
- TSKgel ODS-100V, 3μm, 4.6mm ID x 15cm
- TSKgel precolumn BSA-ODS-100V, 2mm ID x 1cm

Instrumentation

Pump: DP-8020 (Tosoh)

Detector: UV-8020, RI-8020 (Tosoh)

Auto-sampler: AS-8020 (Tosoh)


Data processing: LC-8020 model 2 (Tosoh)

Chemicals and Reagents

All saccharides and reagents were purchased from Kishida Chemicals (Osaka). The pyridylaminated sugars were from Takara Bio Inc. All organic solvents were of HPLC grade from Kishida Chemicals. Human serum was purchased from Sigma.



Figure 1: Separation of Sugars on a TSKgel NH₂-100, 3 μ m Column



TSKgel NH₂-100, 3 μ m, 4.6mm ID x 15cm

Eluent: A) H₂O/Acetone (25/75)
B) 100mmol/L TEA-FA (pH 10.0) /Acetone (25/75)

Flow rate: 1.0mL/min

Injection vol.: 10 μ L

Temperature: 50°C

Detection: RI

Samples: 1. fructose 2. sorbitol 3. glucose

4. sucrose 5. maltose 6. lactose

Table 1: Effect of Organic Solvent on Elution Profiles of Various Sugars

No.	Sample	TEA-FA/ACN			TEA-FA/Acetone			TEA-FA/EtOH		
		Rt (min)	Area	TP	Rt (min)	Area	TP	Rt (min)	Area	TP
1	Rhamnose	3.74	516	9,773	2.98	508	4,662	3.10	413	1,051
2	Ribose	3.84	319	13,951	3.10	367	14,511	3.38	117	7,181
3	Erythritol	4.01	784	15,498	3.13	732	15,495	3.04	756	11,041
4	Fucose	4.07	674	2,640	3.36	658	1,120	3.64	449	274
5	Xylose	4.42	493	6,463	3.43	569	7,349	3.76	274	5,189
6	Arabinose	4.59	381	6,520	3.62	424	4,026	4.09	79	1,213
7	Xylitol	4.76	831	15,677	3.54	794	15,409	3.58	808	10,820
8	Fructose	5.01	791	14,774	3.75	786	11,165	4.28	794	4,925
9	Sorbose	5.11	858	15,455	3.76	828	14,483	4.16	836	9,043
10	Sorbitol	5.69	830	15,269	3.99	827	15,185	4.27	831	10,113
11	Mannose	5.51	462	10,727	3.94	518	5,517	4.66	247	939
12	Mannitol	5.84	825	15,400	4.05	802	15,321	4.41	797	9,749
13	Galactose	6.02	470	4,125	4.35	596	1,924	5.16	176	443
14	Glucose	5.95	852	6,645	4.25	862	4,508	4.96	759	1,152
15	Sucrose	8.47	887	15,707	5.18	887	14,543	5.96	884	7,206
16	Inositol	8.79	938	12,627	6.35	972	13,576	8.30	909	3,089
17	Cellobiose	9.67	666	8,656	6.89	783	4,786	7.63	572	1,186
18	Maltose	9.78	697	7,448	5.91	788	4,505	7.50	700	1,217
19	Lactose	10.28	635	9,110	6.46	703	5,362	8.67	542	4,512
20	Trehalose	10.53	752	15,164	6.36	797	14,181	8.51	830	742
21	Gentiobiose	11.65	561	7,449	7.04	707	3,923	10.98	297	11,041
22	Glucosamine	6.32	675	8,546	4.14	703	1,950	3.83	610	213
23	N-Acetylglucosamine	5.04	898	1,759	3.70	936	792	3.55	894	145
24	Galactosamine	6.53	621	4,614	4.14	594	926	3.81	526	167



Table 1: Chromatographic Conditions

TSKgel NH₂-100, 3 μ m, 4.6mm ID x 15cm

Eluent: A,B: 100mmol/L TEA-Formic acid (pH10.0) /solvent (25/75)

C: 100mmol/L TEA-Formic acid (pH10.0) /solvent (5/95)

Flow rate: 1.0mL/min

Injection vol.: 10 μ L


Temperature: 40°C

Detection: RI

Samples: saccharides



Figure 2: Comparison of Separation of PA Sugars on RPC & HILIC Columns



A) TSKgel NH₂-100, 3μm, 4.6mm ID x 15cm

B) TSKgel ODS-100V, 3μm, 4.6mm ID x 15cm

Eluent: A) A: 200mmol/L TEA-Acetic acid (pH6.5)/ACN (30/70)
B: 500mmol/L TEA-Acetic acid (pH6.5) /ACN (60/40)

B) A: 50mmol/L HCO₂NH₄/ACN (98/2),
B: 50mmol/L HCO₂NH₄/ACN (90/10),

Gradient: 0min (0% B) 30min (100% B) 45min (100% B)

Flow rate: 1.0mL/min

Injection vol.: 10μL

Temperature: 40°C

Detection: Fs (Ex. 315nm, Em. 380nm)

Samples: pyridylaminated (PA) sugars



TOSOH

Figure 2: Pyridylaminated Sugars





Figure 2: Pyridylaminated Sugars Continued







Figure 3: Separation of Vitamins on TSK-GEL HILIC & ODS Columns



A) TSKgel NH₂-100, 3 μ m, 4.6 mm ID x 15 cm

B) TSKgel ODS-100V, 3 μ m, 4.6 mm ID x 15 cm

Eluent: A) 25 mmol/L Phosphate buffer (pH 2.5)/AcCN=30/70

B) 25 mmol/L Phosphate buffer (pH 2.5)/AcCN=88/12

Flow rate: 1.0mL/min

Injection vol.: 5 μ L

Temperature: 40°C

Detection: UV@254nm


Samples: 1.nicotinamide 2.vitamin B2 3.pyridoxine

4.nicotinic acid 5.vitamin C 6.vitamin B1

7.vitamin B12



Figure 4: Comparison of Elution Pattern of Theophylline and its Metabolites on TSK-GEL HILIC & ODS Columns



A) TSKgel NH₂-100, 3 μ m, 2.0mm ID x 15cm

B) TSKgel ODS-100V, 3 μ m, 2.0mm ID x 15cm

Eluent: A) A: 0.1mol/L Triethylamine-Formic acid (pH 10.0)/ACN (5/95)
B: 0.1mmol/L Triethylamine-Formic acid (pH 10.0)/ACN (50/50)

B) A: H₂O/ACN (98/2)+0.1% Formic acid
B: H₂O/ACN (50/50)+0.1% Formic acid

Flow rate: 0.25mL/min

Injection vol.: 10 μ L


Temperature: 40°C

Detection: UV@254nm

Samples: 1. theophylline 2. 3-methylxanthine
3. 1,3-dimethyluric acid 4. 1-methyluric acid
(concentration: 10 μ g/mL for each)



Figure 5: Separation of Theophylline and its Metabolites in Serum with On-line Deproteinization



TSKgel NH₂-100, 3 μ m, 4.6mm ID x 15cm

TSKgel precolumn BSA-ODS-100V, 2.0mm ID x 1cm

A: Standard B: Serum spiked with the standard samples

Eluent: Pretreatment: 0.2mol/L HCO₂NH₄ (pH 3.6) during 0-0.3min

A: ACN B: H₂O/ACN=15/85 C: 0.2mol/L HCO₂NH₄ (pH3.6)/ACN=30/70 (Stepwise gradient: 0.3-2.0min A, 2.0-8.0min B, 8.0-20min C)

Flow rate: 1.0mL/min Detection: UV@254nm Temperature: 40°C Injection vol.: 5 μ L


Samples: 1. caffeine (Cf) 2. theobromine (Tb) 3. theophylline 4. 3-methylxanthine (3-MX) 5. 1-methylxanthine (1-MX)

6. 1,3-dimethyluric acid (DMU) 7. 1-methyluric acid (1-MU)

50 μ g each



Figure 6: Separation of Theophylline and its Metabolites in Serum with Off-line Deproteinization



TSKgel NH₂-100, 3 μ m, 2.0mm ID x 15cm

A: Supernatant of serum deproteinized with 10-fold ACN

B: Supernatant of spiked serum deproteinized with 10-fold ACN

Eluent: A: 0.1mol/L Triethylamine-Formic acid (pH10.0)/ACN (5/95)

B: 0.1mol/L Triethylamine-Formic acid (pH10.0)/ACN (50/50)

Gradient: 0 min-2min (0% B) 2-30 min (80% B) 30-32min (0% B)

Flow rate: 0.25mL/min

Injection vol.: 10 μ L

Temperature: 40°C

Detection: UV@254 nm

Samples: 1. theophylline 2. 3-methylxanthine 3. 1,3-dimethyluric acid

4. 1-methyluric acid

(concentration: 10 μ g/mL for each)



Conclusions

- Organic solvent effect on the retention behavior of sugars using a TSKgel NH₂-100, 3µm HILIC column was investigated. Ethanol and acetone were utilized for an organic modifier instead of acetonitrile. Lower concentration of ethanol was needed to achieve similar retention times because the polarity was higher than that of the other solvents.
- Although the elution order of pyridylaminated (PA) sugars on the TSKgel NH₂-100, 3µm HILIC column was different from that on an ODS column, two acidic PA sugars were completely separated on the TSKgel NH₂-100 column due to an additional ionic interaction. Many water soluble vitamins eluted earlier from the ODS column with poor resolution, whereas excellent separations on the TSKgel NH₂-100 HILIC column were obtained due to stronger retention for hydrophilic compounds.
- On-line drug monitoring was investigated. Theophylline and its metabolites were separated by using a column switching technique. However, hydrophobic compounds which cannot be strongly retained, eluted from the column with broader peak shapes at high concentrations of ACN. Trace amounts of samples were carried over to the next analysis. Further optimization of the operating conditions would be necessary.
- In the case of the off-line deproteinization, the spiked samples could be separated successfully. A TSKgel NH₂-100, 3µm HILIC column can be useful for analyzing hydrophilic compounds in blood and is applicable to the separation of other hydrophilic drugs.