



Reduce Downstream Processing Costs for mAbs by Switching to a Tosoh 2-Step Platform

Introduction

Downstream processing is responsible for up to 80 percent of the entire production costs of biopharmaceuticals. Given the current drive to reduce the cost of manufacturing for biological therapeutics, streamlining downstream processing is a necessity for chromatographers and process engineers.

In this study, we showcase the benefits of a Tosoh 2-step process for the purification of monoclonal antibodies (mAb) in comparison to the standard industrial process. Combining high performance protein A capturing and a single polishing step on a salt-tolerant anion exchange resin (AEX), we could reduce the downstream costs by 45% and increase the production output by 58%.

Materials and Methods

Resins and Pre-packed Columns

TOYOPEARL® AF-rProtein A HC-650F is a high capacity protein A resin for the purification of mAbs. This resin exhibits dynamic binding capacities (DBC) of 70 g/L at 5 minutes of residence time.

TOYOPEARL NH₂-750F, a salt-tolerant anion exchange resin, is based on the TOYOPEARL methacrylate backbone and is functionalized with primary amine groups. TOYOPEARL NH₂-750F resin is ideal for the intermediate purification of mAbs and other proteins. Impurities, such as DNA, viruses, host cell proteins, and endotoxins, are removed. Furthermore, because of the relatively low pK_a value (between 7 and 9), TOYOPEARL NH₂-750F is also able to remove mAb aggregates efficiently. With this unique feature, both polishing steps that are usually necessary to remove all impurities are combined in one flow-through polishing step.

All experiments were performed on SkillPak™ 5 mL pre-packed columns. The SkillPak columns are designed for fast method development or resin screening. These columns guarantee optimal performance and can be operated with commonly used low or medium pressure liquid chromatography systems. They are reproducibly packed and take into account the varying compressibility of each resin, providing an accurate representation of conditions found in full-scale columns.

Purification Protocols

Capture – TOYOPEARL AF-rProtein A HC-650F

TOYOPEARL AF-rProtein A HC-650F was equilibrated with 100 mmol/L sodium phosphate, pH 7.0, and loaded with 20 CV (100 mL) of clarified cell culture fluid with 2 mg/mL Adalimumab. The washing step was performed with 100 mmol/L sodium acetate, pH 7.0, for 10 CV (50 mL). Elution was carried out with 100 mmol/L sodium acetate, pH 3.0. The cleaning of the column was performed with 200 mmol/L NaOH. After cleaning, the column was re-equilibrated with sodium phosphate buffer. The flow rate in the steps of equilibration, washing, elution, and re-equilibration was 204 cm/hr, 150 cm/hr for load, and 180 cm/hr for CIP.

Virus Inactivation

The eluate from the capture step was held at pH 3.0 for 1 hour before being adjusted back to pH 8.0 and the desired conductivity for the polishing experiments with 1 mol/L Tris. The aggregate content after the incubation at pH 3.0 was 1.05%.

Polish – TOYOPEARL NH₂-750F

The protein A purified antibody was diluted to a concentration of 1 mg/mL and adjusted to the conductivity of 20, 22.5 or 25 mS/cm and 20 mmol/L Tris-HCl, pH 8. The column was equilibrated with different conductivities (20, 22.5, and 25 mS/cm – adjusted with NaCl) of the equilibration buffer. The sample was loaded for 40 CV (200 mL). Afterward, a washing step with equilibration buffer for 5 CV (25 mL) was performed. The cleaning of the column was carried out with 500 mmol/L NaOH. The flow rate during the entire process was 300 cm/hr.

Analytical SEC

The load after the 1 hour hold at pH 3.0 and the flow-through of TOYOPEARL NH₂-750F was analyzed by size exclusion chromatography (SEC) using a TSKgel® UP-SW3000 column with the following conditions:

| | |
|-----------------|---|
| Mobile phase: | 100 mmol/L phosphate buffer + 100 mmol/L sodium sulfate + 0.05% NaN ₃ , pH 6.7 |
| Flow rate: | 0.35 mL/min |
| Temperature: | 25°C |
| Detection: | UV @ 280 nm |
| Injection vol.: | 10 µL |

Results and Discussions

Process Development

In previous work, we optimized the capture of the antibody on the TOYOPEARL AF-rProtein A HC-650F. After the capture step, we introduced a virus inactivation step by holding the eluate for 1 hour at low pH.

The first experiments on the salt-tolerant anion exchanger TOYOPEARL NH₂-750F were performed in bind-and-elute mode with a linear gradient at pH 8.0 to determine the necessary conductivity for the flow-through experiments. (Figure 1).

As the monomer eluted first, the next step consisted in adapting the method to achieve flow-through purification. For the flow-through experiments, we used conductivities between 20 and 25 mS/cm, as the monomer is eluting under these conditions, whereas the aggregates remain bonded to the resin. The purities and yields at three different conductivities are listed in Table 1.

Figure 1. Bind- and elute chromatogram of a mAb purification on SkillPak TOYOPEARL NH₂-750F 5 mL column





Table 1. Aggregate content and yield of the flow-through process on TOYOPEARL NH₂-750F at different conductivities

| Conductivity | Aggregate | Yield |
|--------------|-----------|--------|
| 20 mS/cm | 0 % | 91.3 % |
| 22.5 mS/cm | 0.09 % | 93.3 % |
| 25 mS/cm | 0.27 % | 94.6 % |

To achieve the desired purity without a subsequent chromatography step, we chose the conductivity of 20 mS/cm for the platform design. The corresponding chromatogram is shown in *Figure 2*.

Figure 2. Flow-through chromatogram of a mAb purification on SkillPak TOYOPEARL NH₂-750F 5 mL column



Purification Platform

Both chromatography steps were now combined in one integrated process, including the intermediate low pH hold. Protein A has a recovery of 98.8%, while AEX has a recovery of 91.3% (20 mS/cm), which results in a total recovery of 90.2%. DNA, HCP and leached protein A were removed to the Limit of Detection of the used assays (see *Table 2*).

TOYOPEARL, TSKgel and Tosoh Bioscience are registered trademarks of Tosoh Corporation. SkillPak is a registered trademark of Tosoh Bioscience LLC in the USA.

Table 2. Critical quality attributes of a Tosoh 2-step purification platform

| Critical Parameter | Aggregate (%) | Yield (%) | DNA (ppm) | Host Cell Protein (ppm) | Leached Protein A (ppm) |
|----------------------|---------------|-----------|-----------|-------------------------|-------------------------|
| Feed | 0.56 | | 5,500 | 550,000 | < 0.5 |
| After protein A | 1.05 | 98.8 | 60 | 19,000 | 5.0 |
| After AEX (20 mS/cm) | 0.00 | 90.2 | <0.2 | < 30 | < 0.5 |

Cost Analysis

We used the BioSolve (Biopharm Service Ltd., UK) bioprocess analysis software to compare the downstream costs of the optimized two-step process with the costs of an industry standard process published by BioPhorum Operations Group. As shown in *Table 3*, the optimized Tosoh process offers 45% lower costs per gram than the standard industrial process. In addition, process times have been reduced by 58% due to the elimination of one chromatography step and higher flow rates on the Tosoh resins.

Table 3. Comparison of the purification processes using BioSolve Software

| | BioPhorum process (3 steps) | Tosoh process (2 steps) |
|-----------------------------------|-----------------------------|-------------------------|
| Batches per year | 213 | 510 |
| Throughput (doses or kg per year) | 717 | 1,793 |
| Total capital (USD M) | 35.8 | 36.8 |
| Calculated Capacity Utilization | 80% | 80% |
| Cost per gram (USD) | 44.49 | 24.45 |
| Capital Charge | 10.38 | 4.13 |
| Materials | 0.91 | 0.37 |
| Consumables | 8.91 | 4.45 |
| Labor | 21.20 | 14.20 |
| Others | 3.09 | 1.29 |
| Cost Breakdown (%) | 100% | 100% |
| Capital Charge | 23% | 17% |
| Materials | 2% | 2% |
| Consumables | 20% | 18% |
| Labor | 48% | 58% |
| Others | 7% | 5% |

Conclusion

TOYOPEARL NH₂-750F is an effective anion exchange resin for the removal of dimer and higher-order aggregates from mAb monomer. This is the ideal resin for streamlining purification processes of monoclonal antibodies in combination with high performance protein A resins.

The Tosoh 2-step antibody purification process presented here using SkillPak pre-packed columns combines high recovery with low process costs and processing time. This platform offers a superior alternative to the industry standard 3-step process. It can easily be scaled up to a pilot plant and eventually to a manufacturing scale for increased productivity and profitability.