



# Octave™ PRO Specification Sheet


## System Overview

The Octave PRO system from Tosoh Bioscience is a cGMP-ready multi-column chromatography (MCC) system for the purification of biopharmaceuticals at the clinical and commercial scale. The Octave PRO is designed for direct scale-up after process development on the Octave BIO, and both systems offer the same method-execution capabilities. The Octave PRO enables ultimate process versatility and agility to address the most challenging downstream bottlenecks.

The entire flow path of the instrument is single-use, easy to replace and gamma irradiated. The flow path features a biocompatible, pneumatically actuated 104-valve array that directs process fluids through 1 to 8 columns accommodating increasingly intensified upstream titers. Fluid flow is controlled by six independent pumps with a maximum flow rate of 150 L/h, which operate based on real-time data feedback from six integrated flow sensors and six pressure transducers. Additionally, the system includes four

sensor towers each containing a dual-channel UV sensor, a conductivity monitor, and a pH sensor housed in a compact single-use cell for fluid stream monitoring. The Octave PRO's three-dimensional design gives it a compact footprint. Pump head fluid connections, HMI keyboard, and touch screen monitor are arranged on the front with column connections to the valve block and outlet line connections to collection vessels on the back.

Tosoh Bioscience's PROController™ software provides automated control of the Octave PRO system using the AVEVA Wonderware Framework for security and data integrity under Good Automatic Manufacturing Practices issue 5 (GAMPV). The application communicates via a supervisory control and data acquisition (SCADA) computer to a programmable logic controller (PLC), which in turn controls the equipment (pumps, valves, sensors, nitrogen supply, power, communication interfaces, etc.). The 21 CFR part 11 compliant software allows for process control in a cGMP environment.



## Technical Specifications - System Architecture

| Description                   | Specifications  |
|-------------------------------|---|
| Dimensions (W / D / H)        | 213 / 115 / 219 cm, (198 cm without alarm tower)              |
| Weight                        | 450 kg  |
| Operating flow rate           | 1 – 150 L/hr (2.5 L/min)                                      |
| Maximum operating pressure    | 6 bar (87 psi)  |
| Valve Pressure                | 8 bar (116 psi)   |
| Gas inlet pressure            | 9 bar (130 psi)   |
| Column Positions/ Connections | 1 to 8 columns, 3/4 in. sanitary clamps                       |
| Valve block                   | 1 single-use assembly containing 104 pneumatic two-way valves |

| Description          | Specifications   |
|----------------------|--|
| Inlet number         | 6  |
| Outlet number        | 6  |
| Pumps                | 6; 4-piston-diaphragm pump, flow rate ranges 1-150 L/h   |
| Flowmeters           | 6; 1 after each pump, flow rate ranges 0.01-8 L/min      |
| Pressure sensors     | 6; 1 after each pump, measuring range 1-10 bar           |
| Outlet sensors       | 4 assignable single-use combination UV/ conductivity/ pH |
| UV sensors           | UV 280 nm & 305 nm                                       |
| Conductivity sensors | CND range 0–850 mS/cm                                    |
| pH sensors           | pH range 0-14  |

## Single-Use Flowkit

| Item              | Description  |
|-------------------|--|
| Inlet Connections | 6 Aseptiquik® G  |
| Outlet Connection | 6 Aseptiquik® G  |
| Pumps             | 6 Quattroflow® EZ-Set Pump Chamber quaternary diaphragm single-use heads |
| Pressure sensors  | 6 single-use gauge tees (TS1, 1 after each pump)                         |
| Flow sensors      | 6 single-use; ultrasonic (1 after each pump)                             |

## Single-Use Materials of Construction

| Component                              | Wetted Materials  |
|--|---|
| Pump head                              | polypropylene (PP), thermoplastic elastomer (TPE), ethylene propylene diene monomer (EPDM)  |
| Pressure adapter                       | PP, TPE   |
| Flow sensor                            | PP  |
| Valve block assembly                   | USP class VI compliant polyetherimide (PEI), polyvinylidene fluoride (PVDF), polyfluoroalkoxy (PFA)   |
| Optical/ conductivity/<br>pH flow cell | quartz (UV-transparent), EPDM, stainless steel 1.4435 (SS 316L), polyphenylsulfone (PPSU)   |
| Fittings                               | PEEK, platinum cured silicone, PP   |
| Tubing                                 | ¼ in. ID polybraided platinum cured silicone, pump inlet is weldable ¼ in. ID Advantaflex biopharmaceutical grade thermoplastic elastomer (TPE) |
| Gamma irradiation dose                 | 25-40 kGy   |

## Electrical Requirements

| Parameter       | North America            | EMEA                              |
|-----------------|--------------------------|-----------------------------------|
| Voltage         | 120V; 1 phase            | 220-240V; 1 phase                 |
| Amps            | 30A                      | 15A                               |
| Frequency       | 60 Hz                    | 50 Hz                             |
| Number of Cords | 1                        | 1                                 |
| Type of Plug    | NEMA L5-30P, with ground | Must meet local code, with ground |

## Ordering Information

### Systems

| P/N     | Description   |
|---------|---|
| 0041200 | Octave PRO US, Octave PRO GMP ready multi-column chromatography skid for US region. |
| 0041201 | Octave PRO EU, Octave PRO GMP ready multi-column chromatography skid for EU region. |
| 0041202 | Octave PRO SU Flowkit, Full single-use flowkit for the Octave PRO.                  |

Tosoh Bioscience is a registered trademark of Tosoh Corporation.  
Octave and PROController are trademarks of Tosoh Bioscience LLC.  
Aseptiquik is a registered trademark of Colder Products Company.  
Quattroflow is a registered trademark of PSG Germany GmbH.