



RESIN INFORMATION SHEET

Product Name TOYOPEARL® Super A
(Functionalised resin for affinity chromatography)

Part Numbers	0023580	TOYOPEARL Super A, 10 mL
	0023581	TOYOPEARL Super A, 25 mL
	0023582	TOYOPEARL Super A, 100 mL
	0023583	TOYOPEARL Super A, 1 L
	0023584	TOYOPEARL Super A, 5 L

Product Description TOYOPEARL chromatographic resins are based on a rigid methacrylic polymer, resulting in high mechanical and chemical stability. Resins are available as non-functionalized “HW” series resins for size exclusion separations, and derivatized with surface chemistries for alternative modes of chromatography such as ion exchange, hydrophobic interaction or affinity separations.

TOYOPEARL Super A resin is a functionalized support resin for affinity chromatography. The resin is prepared by multipoint attachment of a recombinant Protein A variant onto TOYOPEARL non-functionalized resin. TOYOPEARL Super A resin is useful for the chromatographic purification of monoclonal antibodies at process scale. The Protein A ligand is not animal derived.

Operating Conditions	Packing pressure	Typically 0.3 MPa
	Shipping solvent	20 % (v/v) ethanol
	Shipping formulation	72 % (v/v) slurry in shipping solvent (*)
	Pressure limiting factor	Depend on column hardware (typically 0.7 MPa)
	Operating linear flowrate	Typically 10 - 600 cm/hour
	Long-term storage conditions	20 % (v/v) ethanol at 2 - 8 °C
	Cleaning-in-place/Sanitization	0.1 - 0.5 mol/L NaOH

Specifications	Particle size distribution (min. 80 % within range)	30 - 60 µm
	Protein adsorption capacity (of human antibody)	Min. 68 g/L - gel
	Bacterial count	Max. 100 CFU/mL
	Endotoxin concentration	Max. 10.0 EU/mL
	Eluable matter	Max. 0.2 % (for dry gel)
	Foreign substance (colored particle)	Max. 6

Additional Information	Appearance	White resin slurry which settles upon standing
	Mean pore diameter (base resin)	100 nm (*)

(*) The value is for reference only, not guaranteed.

Lot-specific data are included in the Certificate of Analysis (COA) shipped with the product. For detailed test procedures please refer to the appropriate Regulatory Support File.