

TSKgel® G3000SW Products

Part Numbers:	08800, TSKgel G3000SW, 8.0 mm ID x 30 cm, Glass	08805, Guard column for TSKgel G3000SW Glass column, 8 mm ID x 4 cm
	05789, TSKgel G3000SW, 7.5 mm ID x 30 cm	05371, Guard column for TSKgel G3000SW columns, 7.5 mm ID x 7.5 cm
	05103, TSKgel G3000SW, 7.5 mm ID x 60 cm	05758, Guard column for for TSKgel G3000SW columns, 21.5 mm ID x 7.5 cm
	06728, TSKgel G3000SW, 21.5 mm ID x 30 cm	06819, TSKgel SW top-off gel, 1g
	05147, TSKgel G3000SW, 21.5 mm ID x 60 cm	

This sheet contains the recommended operating conditions and the specifications for TSKgel G3000SW columns. Installation instructions and column care information are described in a separate Instruction Manual.

A. OPERATING CONDITIONS	
1. Shipping Solvent:	0.05% NaN ₃ and 0.1 mol/L Na ₂ SO ₄ in 0.1 mol/L phosphate buffer, pH 6.7
2. Max. Flow Rate:	0.8 mL/min (8.0 mm ID Glass) 1.2 mL/min (7.5 mm ID) 8.0 mL/min (21.5 mm ID)
	When a buffer with high viscosity is used, the maximum flow rate may have to be reduced so as not to exceed the maximum pressure drop. When changing solvents, use a flow rate equal to 25% of the maximum flow rate.
3. Standard Flow Rate:	0.4 - 0.8 mL/min (8.0 mm ID Glass) 0.5 - 1.0 mL/min (7.5 mm ID) 3.0 - 6.0 mL/min (21.5 mm ID)
4. Max. Pressure:	1.5 MPa (21.5 mm ID x 30 cm) 2.0 MPa (8.0 mm ID Glass x 30 cm) 2.5 MPa (7.5 mm ID x 30 cm) 3.0 MPa (21.5 mm ID x 60 cm) 5.0 MPa (7.5 mm ID x 60 cm)
5. pH Range:	2.5 - 7.5
6. Salt Conc.:	≤ 0.5 mol/L
7. Organic Conc.:	0 - 100% for aqueous soluble organic solvents. Make gradual solvent changes using a shallow gradient at low flow rate.
8. Temperature:	10 – 30 °C, Reduce flow rate when operating below 10 °C.
9. Cleaning Solvents:	1. Turn the column in reverse flow direction and run at half the maximum flow rate. 2. Clean with 5 column volumes (CV) of 1 mol/L sodium chloride, pH 7.0 3. Clean with 5CV of ultra-pure water. 4. Clean with 5CV of 20% acetonitrile. 5. Clean with 5CV of ultra-pure water. 6. Turn column in normal flow direction and equilibrate in mobile phase for at least 45 minutes
10. Storage:	Store the column in mobile phase containing 0.05% NaN ₃ or 20% ethanol when it will not be used the next day. For overnight storage flush the column with mobile phase at low flow rate. Prevent air from entering the column!
11. Column Protection:	The use of guard columns is recommended to prolong the life of the analytical column. Guard column life depends greatly on sample cleanliness. As a general rule, guard columns should be replaced after every 30-40 sample injections, when the peaks become excessively wide, or when the peaks show splitting.
12. TSKgel SW top-off gel:	Occasionally, due to accident, sample, mobile phase or operational variables, a depression can develop at the column or guard column inlet. Use TSKgel SW top-off gel for filling in such voids.
B. SPECIFICATIONS	
The performance of TSKgel G3000SW columns is tested under the conditions described in the Data Sheet. All columns have passed the following quality control specifications:	
1. Number of Theoretical Plates (N):	≥ 10,000 (other 30 cm columns) ≥ 20,000 (60 cm columns)
2. Asymmetry Factor (AF):	0.7 - 1.6

DS1006 Revised 06APRIL2014

Tosoh Bioscience LLC
156 Keystone Drive
Montgomeryville, PA 18936
Phone: (215) 283-5000
Orders & Technical Service: (800) 366-4875
FAX: (215) 283-5035

Tosoh Bioscience GmbH
Zettachring 6
D-70567 Stuttgart
Germany
Phone: 0711/13257 0
FAX: 0711/13257 89