

## TSKgel® HT and HT2 (High Temperature) Products

<b>Part Numbers:</b>	07112, TSKgel GMHxl-HT, 13 µm, 7.8 mm ID x 30 cm	22889, GMHhr-H(S) HT2, 13 µm, 7.8 mm ID x 30 cm
	18393, TSKgel GMHhr-H(S) HT, 13 µm, 7.8 mm ID x 30 cm	22887, GMHhr-H(30) HT2, 30 µm, 7.8 mm ID x 30 cm
	18391, TSKgel GMHhr-H(30) HT, 30 µm, 7.8 mm ID x 30 cm	22888, GMHhr-H(20) HT2, 20 µm, 7.8 mm ID x 30 cm
	18392, TSKgel GMHhr-H(20) HT, 20 µm, 7.8 mm ID x 30 cm	22890, G2000Hhr HT2, 20 µm, 7.8 mm ID x 30 cm
	18395, TSKgel G2000Hhr HT, 20 µm, 7.8 mm ID x 30 cm	22892, Guard column for 7.8 mm ID TSKgel GMHhr-H(S) HT2, 13 µm, 7.5 mm ID x 7.5 cm
	18420, TSKgel GMHhr-H HT, 5 µm, 7.8 mm ID x 30 cm	22891, Guard column for 7.8 mm ID TSKgel GMHhr-H(30) HT2 and GMHhr-H(20) HT2 columns, 30 µm, 7.5 mm ID x 7.5 cm
	18396, Guard column for 7.8 mm ID TSKgel GMHhr-H(30) HT and GMHhr-H(20) HT columns, 30 µm, 7.5 mm ID x 7.5 cm	
	18397, Guard Column for 7.8 mm ID TSKgel GMHxl-HT & GMHhr-H(S) HT columns, 30 µm, 7.5 mm ID x 7.5 cm	

This sheet contains the recommended operating conditions and the specifications for TSKgel HT and HT2 (High Temperature) columns and guard columns. Installation instructions and column care information are described in a separate Instruction Manual.

<b>A. OPERATING CONDITIONS</b>	
1. Shipping Solvent:	o-dichlorobenzene (ODCB)
2. Standard Flow Rate:	0.5 - 1.0 mL/min
3. Max Flow Rate:	2.5 mL/min at ≤140°C, 1.0 mL/min at >140°C
4. Max. Pressure:	1.5 MPa (GMHxl-HT, GMHhr-H(30) HT, GMHhr-H(30) HT2, GMHhr-H(20) HT, GMHhr-H(20) HT2, G2000Hhr HT, G2000Hhr HT2) 2.0 MPa (GMHhr-H(S) HT, GMHhr-H(S) HT2) 3.5 MPa (GMHhr-H HT)
5. Multiple Columns:	Columns of the same or different pore size are often connected in series to improve resolution and/or to expand the linear portion of the calibration curve. Connect the columns in order of decreasing pore size to avoid overloading from the high MW components. Connect analytical columns using short pieces of 1/16" x 0.01" ID stainless steel tubing.
6. Solvents	TSKgel HT and HT2 columns are only available packed in o-dichlorobenzene (ODCB).
7. Temperature:	TSKgel HT: up to 140°C TSKgel HT2: up to 220°C  (Heat the column gradually to prevent damage. For example, it is recommended to ramp up the temperature from ambient to 140 °C at less than 50 °C/hour or to 220 °C at less than 25 °C/hour.)
8. Sample Size:	0.001 - 0.5 mg (analytical)
9. Storage:	TSKgel HT and HT2 columns can be left overnight in solvent in the LC system. When it will not be used for longer periods of time, remove the column from the equipment, seal the ends with the provided protective screws, and store it at laboratory temperature. At all times, prevent air from entering the column!
10. Column Protection:	The use of guard columns is recommended to prolong the life of the analytical column. Guard columns are not suited for analysis. Their main function is to protect the analytical column; they do not improve resolution when connected to the main column. They are also not a substitute for filtering the mobile phase and the sample. A guard column does reduce pump pulsation, and further protects the main column by collecting highly adsorptive components and insoluble substances. Guard column life depends greatly on sample cleanliness. As a general rule, guard columns should be replaced when the peaks become excessively wide, or when the peaks show splitting.

DS1240 29JULY2016

[www.tosohbioscience.com](http://www.tosohbioscience.com)

**Tosoh Bioscience LLC**  
3604 Horizon Drive, Suite 100  
King of Prussia, PA 19406  
Phone: (484) 805-1219  
Orders & Technical Service: (800) 366-4875  
FAX: (610) 272-3028

**Tosoh Bioscience GmbH**  
Zettachring 6  
D-70567 Stuttgart  
Germany  
Phone: 0711/13257 0  
FAX: 0711/13257 89

<b>B. SPECIFICATIONS</b>	
The performance of TSKgel HT and HT2 (High Temperature) columns is tested under the conditions described in the data sheet. All columns have passed the following quality control specifications:	
1. Number of Theoretical Plates (N):	> 5,500 (GMHxl-HT) > 8,000 (GMHhr-H(S) HT, GMHhr-H(S) HT2) > 4,000 (GMHhr-H(30) HT, GMHhr-H(30) HT2) > 6,000 (GMHhr-H(20) HT, GMHhr-H(20) HT2, G2000Hhr HT, G2000Hhr HT2) > 16,000 (GMHhr-H HT)
2. Asymmetry Factor (AF):	0.7 – 1.6 (all columns)

DS1240 29JULY2016

[www.tosohbioscience.com](http://www.tosohbioscience.com)

**Tosoh Bioscience LLC**  
 3604 Horizon Drive, Suite 100  
 King of Prussia, PA 19406  
 Phone: (484) 805-1219  
 Orders & Technical Service: (800) 366-4875  
 FAX: (610) 272-3028

**Tosoh Bioscience GmbH**  
 Zettachring 6  
 D-70567 Stuttgart  
 Germany  
 Phone: 0711/13257 0  
 FAX: 0711/13257 89